
ROSCOE FROG

COMPRESSOR NEBULIZER

Manual

MODEL: 50000
FROG-TRU
FROG-TRUWB
FROG-TRUWB-PB

Read this manual before operating the nebulizer.
Save these instructions for future reference.

TABLE OF CONTENTS

	Page
1. Introduction	2
2. Illustration of Parts	3
3. Important Safety Precautions And Warnings	4 – 7
4. Getting Started	8 – 10
5. Operating the Compressor Nebulizer	11 – 12
6. Cleaning	12 – 13
7. Maintenance	14
8. Specifications	15
9. Limited Warranty	16
10. Service and Repair	16

INTRODUCTION

Thank you for purchasing the Roscoe Medical Frog Nebulizer Compressor. It is a compact medical device designed to efficiently deliver physician prescribed medication to the bronchial lung passages and is intended for home health care use. With proper care, it will provide you with many years of reliable treatment.

This product is developed for the successful treatment of asthma, allergies and other respiratory disorders. It creates a stream of air that travels through a clear tube to the nebulizer. When air enters the nebulizer, it will convert the prescribed medication into an aerosol mist for easy inhalation.

Your nebulizer compressor should be used under the supervision of a licensed physician and/or a respiratory therapist. We encourage you to thoroughly read this manual to learn about the features of this product. Any use of this product other than its intended use should always be avoided.

All Rights Reserved
42-50000_04 ©2018


ILLUSTRATION OF PARTS

NEBULIZER COMPRESSOR UNIT

1. Compressor
2. Cap/Air Filter
3. Air Outlet
4. Power Switch

Figure 1

NEBULIZER ASSEMBLY AND ACCESSORIES ILLUSTRATION

1. Tubing
2. Nebulizer Cup
3. Nebulizer Insert
4. Nebulizer Cap
5. Mouthpiece
6. Child Mask


Figure 2

SAFETY SYMBOLS USED IN THIS MANUAL	
 DANGER	Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.
 WARNING	Indicates a potentially hazardous situation which, if not avoided, could result in serious injury and equipment damage.
 CAUTION	Indicates a potentially hazardous situation which, if not avoided, could result in minor or moderate injury to the user or patient or damage to the device or other property.

IMPORTANT SAFETY PRECAUTIONS AND WARNINGS

It is important that you read all warnings in user's manual before use. Use the product only for its intended use as described in the manual. Improper use can result in serious or fatal illness / injury, improper treatment or property damage. **DO NOT** attempt to operate this product without reading all instructions carefully. Save this manual for future reference.

 WARNING (USA only)	U.S. Federal Law restricts these devices to sale by, or on the order of, a physician or licensed practitioner. This device should only be used under the continued supervision of a physician or licensed practitioner.
 DANGER	THIS MEDICAL DEVICE IS NOT INTENDED TO BE USED AS A LIFE SUSTAINING OR LIFE SUPPORTING DEVICE.

This device has been thoroughly tested and inspected to assure proper performance and operation.

Specifications put forth in this manual were in effect at the time of publication. However, to ensure continual improvement measures, changes to these specifications may be made at any time without obligation on the part of manufacturer.

Upon receiving your nebulizer/compressor, inspect the packaging for signs of damage in transit. Inspect the unit thoroughly. If damage is noted, external or internal, **DO NOT** use, and report to your Homecare Dealer immediately.

When using electrical products, especially when children are present, basic safety precautions should always be followed. Read all instructions before using.

 **WARNING**

To reduce the risk of serious or fatal injury from electrocution, fire, burns and to reduce the risk of damage or malfunction of the unit:

- Before use, check the serial number label on the bottom of the compressor to ensure that the voltage and current indicated on the unit correspond with the voltage and current available. **DO NOT** overload wall outlets or use with power strips or extension cords.
- **DO NOT** use this product with a DC to AC converter, or with any AC voltage and/or frequency other than specified
- Before use, check tubing and accessories for proper assembly. All parts should be firmly in place. Use of improperly assembled tubing and accessories could diminish or prevent adequate delivery or medication and the effectiveness of the treatment.
- Tubing and accessories are for single patient use only. **DO NOT** share your tubing and accessories with other patients.
- **NEVER** use tubing with a clogged nozzle. If the nozzle is clogged, the aerosol mist will be reduced, altered, or stopped, which will diminish or eliminate the effectiveness of the treatment. If clogging occurs, replace with the same model.

DISCONTINUE OPERATION IMMEDIATELY IF ANY OF THE MOTOR HOUSING PARTS BECOME DETACHED OR BROKEN, EXPOSING THE MOTOR OR ANY OTHER ELECTRICAL COMPONENTS.

- **DO NOT** open or disassemble the compressor. Refer all servicing to a Compass Health Brands authorized dealer.
- **ALWAYS** unplug immediately after each use and keep unplugged when not in use. **DO NOT** run this product unattended.
- **NEVER** operate this product if it has a damaged cord or plug, is not working properly or any abnormalities occur, has been dropped in any way or into water, or has been exposed to any liquids inside the case. Return to a Compass Health Brands authorized dealer for repair or replacement.
- **DO NOT** use the nebulizer while operating a vehicle or while bathing.
- Keep all electrical cords away from heated surfaces. **DO NOT** operate or store unit under direct sunlight, high temperature or humidity.
- The compressor motor has a thermal protector, which will shut off the unit before it becomes overheated. When the thermal protector shuts the unit off, please:
 - Switch off the unit and unplug the unit from the electrical outlet.
 - Wait 40 min. for motor to cool down before starting another treatment.
 - Make sure the air openings are not obstructed and air filter is clean.
- **DO NOT** operate where oxygen is being administered in a closed environment such as an oxygen tent.

 **WARNING (cont.)**

- To avoid EMI/RFI:
 - a. **DO NOT** turn on or use hand-held personal electronic communication devices such as cellular phones, walkie-talkies and CB radios while the compressor is turned on.
 - b. **DO NOT** operate the compressor close to any nearby transmitters (radio, television, microwave, etc.).
 - c. Be aware that adding accessories or components or modifying the compressor may make it more susceptible to EMI/RFI.
- **NEVER** block air openings of the main unit or place it where the air openings may be obstructed.
- **NEVER** drop or insert any object into any openings.
- **DO NOT** wrap cord around the appliance.
- **NEVER** allow the power supply cord to be pulled, jerked, strained, twisted, or severely bent, especially at the plug connections. Damage will occur at the high flex point of entry into the appliance, causing it to rupture and short. Handle the cord carefully for longer life.
- Close supervision is necessary when this product is used by or near children and physically or mentally impaired persons. **NEVER** use when sleeping or drowsy.
- This product contains small parts that may present a choking hazard to small children. The tubing also presents a strangulation hazard. Keep the unit out of reach of small children.

WARNING: ALWAYS USE CLOSE ADULT SUPERVISION WHEN ADMINISTERING TREATMENT TO A CHILD. DO NOT LEAVE CHILD UNATTENDED DURING TREATMENT.

- Children under the age of 3 or any patient who is unable to use a mouthpiece properly under supervision should use a mask.
- **DO NOT** place or store the product where it can fall or be knocked into a bathtub, sink, water, or other liquid. **DO NOT** immerse or drop into water or other liquid.
- **DO NOT** reach for a product that has fallen into water or other liquid. Unplug immediately.
- **ALWAYS** disconnect the unit from the electrical outlet before cleaning.
- To reduce the risk of increased bacterial growth, infection, illness or injury from contamination, thoroughly clean and dry all parts of tubing and accessories and thoroughly dry any moisture or condensation in the tubing at the end of every treatment, following the instructions on page 12.
- **NEVER** spray liquids onto the compressor housing. Fluid could cause damage to the electrical compressor parts and could lead to a malfunction. If fluids enter the unit, return the product to an authorized medical equipment service center for examination and repair.

 **WARNING (cont.)**

- Disposable nebulizers are **NOT** dishwasher safe.
- Use a clean, damp cloth to wipe the exterior of the compressor housing.
- Cleaning the nebulizer after each treatment is not enough. The tubing and accessories **MUST** also be regularly sterilized and disinfected between treatments to prevent serious or fatal illness caused by contamination.

 **CAUTION**

OPERATING CAUTIONS:

- Connect this product to an appropriate voltage outlet for your model. **DO NOT** use with power strips or extension cords.
- **DO NOT** use this product with a DC to AC converter, or with any AC voltage and/or frequency other than specified.
- **DO NOT** run this product unattended.
- The compressor motor has a thermal protector, which will shut off the unit before it becomes overheated. When the thermal protector shuts the unit off, please:
 - a. Switch off the unit and unplug the unit from the electrical outlet.
 - b. Wait 40 min. for motor to cool down before starting another treatment.
 - c. Make sure the air openings are not obstructed and air filter is clean.
- **NEVER** operate if this unit has a damaged cord or plug, if it has been dropped in any way or into water, if it does not work properly. Return it to a service center for repair.
- If any abnormality occurs, discontinue use immediately until the unit has been examined and repaired.
- **ALWAYS** unplug the product immediately after use.
- **NEVER** allow the power supply cord to be pulled, jerked, strained, twisted or severely bent, especially at the plug connections. Damage will occur at the high flex point of entry into the appliance, causing it to rupture and short. Handle the cord carefully for longer life.
- **NEVER** block the air openings of the main unit or place it where the air openings may be obstructed.

STORAGE CAUTIONS:

- **DO NOT** store the unit under direct sunlight, high temperature or humidity.
- Keep the unit out of reach of small children.
- Keep the cord away from heated surfaces. Do not wrap the cord around the appliance.
- **ALWAYS** keep the unit unplugged while not in use.

CAUTION (cont.)

CLEANING CAUTIONS:

- **DO NOT** immerse the unit in water. It may damage the unit.
- Disconnect the unit from the electrical outlet before cleaning.
- Clean all necessary parts after each use as instructed in this guidebook.
- Disposable nebulizers are **NOT** dishwasher safe

CLEANING THE COMPRESSOR HOUSING:

- Use a clean, damp cloth to wipe the exterior of the compressor housing.
- **NEVER** spray liquids onto the compressor housing. Fluid could cause damage to the electrical and compressor parts, and could lead to a malfunction. In the event that fluids enter the unit, return the product to an authorized medical equipment center for examination and repair.

GETTING STARTED

1. Check to ensure the following components of the nebulizer are enclosed:
 - Compressor unit
 - Tubing and accessories (see page 3 for illustration of parts)
2. It is recommended to keep extra tubing, accessories and air filters on hand at all times in case original parts become lost, damaged, or soiled.

CAUTION


If you did not receive the components listed in this manual, please contact the dealer where you purchased your nebulizer.

SETTING UP THE NEBULIZER COMPRESSOR AND NEBULIZER ACCESSORIES

CAUTION

Tubing and accessories are **NOT** sterile from package. Prior to initial operation, the nebulizer should be thoroughly cleaned.

1. Wash hands with soap and warm water and dry completely with a clean towel.
2. Place the compressor on a clean, dust-free, stable and flat surface. Check the filter monthly or every 10 – 12 treatments and replace if dirty, clogged, or wet to ensure proper operation.
3. Make sure accessories have been cleaned and/or disinfected before first use and after last daily treatment. Refer to the Cleaning and Disinfecting Tubing and Accessories section on page 12.
4. Remove the nebulizer cap by twisting counter-clockwise until it stops and lift up (Figure 3-1). Be sure that the cone is put inside the bottom section
5. Fill nebulizer cup using eyedropper or syringe with prescribed medication. Do not fill more than 5 ml per treatment (Figure 3-2).
6. Close the nebulizer cap and turn clockwise until it is tightened (Figure 3-3).

7. If using the mask (optional), press it directly into the nebulizer cap. (Figure 4)
8. If using the mouthpiece, attach it firmly into the nebulizer cap. The cap is designed to fit tightly (Figure 5).


Figure 4


Figure 5

9. Attach one end of the tubing firmly onto to the air outlet on the back of the compressor unit. It is important that the tubing be pushed all the way onto the air outlet. The tubing is designed to fit tightly (Figure 6-1).
10. Firmly press the other end of the tubing onto the bottom of the nebulizer cup (Figure 6-2).


Figure 6-1


Figure 6-2

11. Plug the power cord into an appropriate electrical outlet. (**DO NOT** use with power strips or extension cords). Make sure at this stage, the power switch is at “OFF” status.
12. Press power switch to “ON” status to begin your prescribed treatment. (**DO NOT** run unattended). When treating children, adult supervision required.
13. When treatment is finished, shut off the unit and unplug it from the electrical outlet.

⚠ WARNING!

- The exterior surfaces of the unit and the air outlet connector may become slightly warm to the touch during use. Warm air can also be felt exiting the vents on the sides of the case of the unit. This is normal.
- Never operate if this unit has a damaged cord or plug, has been dropped in any way, immersed into water or if it does not work properly. Return it to a service center for repair.
- If any abnormality occurs, discontinue use immediately until the unit has been examined and repaired.
- Always unplug the product immediately after use.
- Never allow the power supply cord to be pulled, jerked, strained, twisted or severely bent, especially at the plug connections. Damage will occur at the high flex point of entry into the appliance, causing it to rupture and short. Handle the cord carefully for longer life.
- Never block the air openings of the main unit or place it where the air openings may be obstructed.
- Safety overload: **DO NOT** exceed 20 minutes of continuous operation.

OPERATING THE COMPRESSOR NEBULIZER

1. Take the power cord and plug into an electrical outlet. (Figure 7)




Figure 7

2. Locate the power switch on the top of the unit. Switch it to the "ON" position as shown below (Fig. 8). Power switch will light up when on.


"ON" Position


Figure 8

3. Check for moisture leakage near the nebulizer connections.
4. There are (3) ways to administer medication with this nebulizer. Consult with your healthcare provider to determine which option is right for you.
 - a. Hold the air intake near the nose (Figure 9-1).
 - b. If using a mask (optional), put the mask on the face (Figure 9-2).
 - c. If using a mouthpiece, hold it in or near the mouth (Figure 9-3).

Nebulizer Unit
Figure 9-1


Mask (optional)
Figure 9-2

Mouthpiece
Figure 9-3

5. Taking deep breaths in, inhale the liquid medication vapor.
6. When liquid medication is depleted, turn the power switch to “OFF” (Fig. 10).

"OFF" Position

Figure 10

⚠ CAUTION

The compressor is designed for intermittent use only.
For optimal performance, operate the unit for 20 minutes on,
then leave off for 40 minutes.

7. Unplug the power cord from the electrical outlet.
8. Detach the tubing and accessories from the compressor unit.
9. Clean nebulizer components before next treatment.

CLEANING

CLEANING THE NEBULIZER COMPONENTS AND ACCESSORIES:

1. Disconnect the tubing from the compressor and from the bottom of the nebulizer cup.
2. If there is any moisture (condensation) in the tubing, let the compressor run with the tubing only for 2 – 3 minutes. The warm air from the compressor will dry out the tubing. Alternatively you can dry the tubing by removing it from the compressor and hanging it up with the “ends down” to allow the moisture or condensation to drain out and evaporate. Dampen a cloth with warm water to wipe the exterior of the tubing.
3. Disassemble tubing end of the accessories.
4. Wash all parts with warm water and liquid dish soap. It is not necessary to wash the tubing.
5. Rinse thoroughly with warm water and shake out water.
6. Air dry or hand dry nebulizer parts with a clean, lint free cloth. Reassemble nebulizer when parts are dry and store. Make sure all components are completely dry before storing.
7. You can also wash the reusable nebulizer (optional accessory) in the top

shelf of a dishwasher. Masks and disposable nebulizers cannot be washed safely in a dishwasher. It is not necessary to wash the tubing. Place the parts in a dishwasher basket and place the basket on the top rack of the dishwasher. Remove and dry parts when the cycle is complete.

8. With proper care, all components of a disposable nebulizer kits should be replaced every 2-4 weeks, and all components of a reusable nebulizer kits should be replaced every 6 months.

 **CAUTION**

Disposable nebulizers are finished. **NOT** dishwasher safe!

CLEANING THE NEBULIZER COMPONENT:

1. Clean the nebulizer components and accessories as described in the cleaning section. Every other treatment day, soak the parts in a solution of one (1) part distilled white vinegar and three (3) parts hot tap water for one (1) hour. **DO NOT SOAK THE TUBING.** You can substitute commercially available respiratory equipment disinfectant for distilled white vinegar by following manufacturer's instructions for mixing.
2. Rinse all parts thoroughly with warm tap water and dry with a clean, lint-free cloth. Discard vinegar solution when disinfecting is complete.

 **WARNING**

Only one patient may use each set of tubing and accessories.
DO NOT share components with other patients.

CLEANING THE COMPRESSOR

1. Wipe daily with a damp cloth.
2. Do not use any powdered cleaners or soap pads, which may damage the finish.

MAINTENANCE

NOTE

The following accessories and replacement parts are recommended

Accessories	Part No.
Kit: Mouthpiece, Disposable Nebulizer Cup, Neb T, 6" Reservoir Tubing & 7' Tubing	NEB-ROSKT
Kit: Adult Mask, Disposable Nebulizer Cup & 7' Tubing	NEB-ADLTM
Kit: Pediatric Mask, Disposable Nebulizer Cup & 7' Tubing	NEB-PEDM
Kit: Mouthpiece, TruNeb Reusable Nebulizer Cup & 7' Tubing	50016
Air Filters (6)	CIF-1128PK6
Nylon Carry Bag	50002

AIR FILTER:

- The air filter on the side panel of the compressor unit is bright white when new. With use, it will turn gray, which means it is time to replace it. It must be replaced immediately if it gets wet or becomes clogged.
- To replace the air filter, ensure the compressor is turned off.

FILTER CHANGE

1. **DO NOT** wash or clean the filter. **ONLY** use filters supplied by your distributor, and **DO NOT** operate without a filter.
2. Check the filter every 15 days. Change the filter when it is visibly dirty OR every 30 days, whichever comes first.
3. Changing procedure:
 - a. Remove the filter cover by twisting counter-clockwise and pulling out.
 - b. Replace the used filter with a new one.
 - c. Replace the filter cover.

! CAUTION

DO NOT attempt to wash and reuse the air filter. **DO NOT** use substitute filter material. Compressor damage may result.


Figure 12

NEBULIZER SPECIFICATIONS

Power	AC 120V, 60Hz
Maximum Wattage	180VA
Power Consumption	1.5A
Medication Capacity	5ml
Particle Size (micron)	0.5 μ m ~ 5 μ m
MMAD	4 μ m
Sound Level	Less than 60 dB at 1 meter
Average Nebulization Rate	\geq 0.2ml/min
Compressor Pressure Range	17.4 to 29 PSI (0.12 to 0.20 MPA)
Operating Pressure Range	8 to 13 PSI (0.05 to 0.07 MPA)
Liter Flow Range	4.0 ~ 8.0 LPM
Operating Temperature Range	10°C to 40°C (50°F to 104°F)
Operating Humidity Range	10 to 90% RH
Storage Temperature Range	-25°C to 60°C (-13°F to 140°F)
Storage Humidity Range	10 to 95 % RH
Dimension (L x W x H)	5.25" x 6" x 7.5" (13.3 cm x 15.2 cm x 19 cm)
Weight	3.3 lbs. (1.5 kg)
Standard Accessories	Nebulizer, Air Tube, Mouthpiece, Child Mask, Air Filters (5pcs)
Optional Accessories	Adult mask, Reusable Nebulizer, Carry Bag
Warranty	5-year limited warranty

LIMITED WARRANTY

Compass Health Brands offers a limited five (5) year warranty on the compressor unit and six (6) months on accessories. We warrant that the compressor unit and its parts will be free from defects in workmanship and materials for a period of five (5) years from the date of original purchase and thirty (30) days on tubing and accessories. This warranty does not extend to failures resulting from accident, misuse, abuse, alteration, use of unauthorized service, parts or cleaning solutions, or failure to comply with the **User's Manual**.

If the warranted equipment should fail during the warranty period, Compass Health Brands, at our option, will repair or replace parts at issue, provided the claim is bona fide. It shall be the responsibility of the purchaser to pack and return the equipment in a manner to avoid shipping damage. The shipment must be prepaid to us. This warranty covers the cost of labor incurred in the removal or replacement of the warranted component parts at us only.

SERVICE AND REPAIR


If you should have any questions about this nebulizer, please contact an authorized Compass Health Brands dealer. Do not open or disassemble the compressor unit at any time. This will void the warranty. If the unit is malfunctioning, discontinue use immediately and return to the original provider for repair or replacement. Please dispose of this nebulizer in accordance with local, municipal, state, and federal law of electronic devices.

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Compass Health Brands Corp.