



Ultra Grip™

Bath suction grab bar Barra de agarre para succión de baño

TO INSTALL:

1. Clean the mounting surface and suction cups with a damp cloth. Rubbing alcohol is recommended to remove residue.
2. Flip levers up
3. Press dual suction cups firmly against a smooth, flat, non porous surface
4. Flip levers down on each end of the suction handle
5. Check that indicator shows fully green (not red)
6. Check that the grip bar is firmly in place before each use, and reinstall if not firm.

NOTE: Suctions cups will stick to surfaces without being locked into place, so always assure that the indicator shows fully green and manually test for adhesion prior to use. Suction handle is for balance only. **DO NOT** place full weight on it or use it to pull yourself up. Compass Health Brands assumes no responsibility for any damage or injury caused by improper installation, assembly, or use of this product.

TO REMOVE:

1. Flip the levers up
2. Lift up using tabs on suction cups to release

WARNINGS: Use caution on wet surfaces. Use to assist with balance only. Grab bars are **NOT** designed to support the total weight of an individual. **DO NOT** use to leverage body weight. **DO NOT** use to pull self out of bathtub. Users with limited physical capabilities should be supervised or assisted in bath and commode areas, even when using the grab bar. Must be mounted on a clean, smooth, flat, non porous surface, and can not be installed over grout lines on tile. Clean the mounting surface and the suction cups before the first use and when necessary. **ALWAYS** check to make sure both indicators are green and test to make sure the bar is securely attached before using. Uninstall and reinstall as needed.