



Deluxe Overbed Table

Model No. PB3037

DEALER: These instructions MUST be given to the user of this product.

SAVE THESE INSTRUCTIONS

NOTE: Check ALL parts for shipping damage. If shipping damage is noted, DO NOT use. Contact Carrier/Dealer for further instruction.

Safety Summary

To ensure the safe and proper installation of the deluxe overbed table, these instructions MUST be followed:

⚠ WARNING

- **DO NOT use this product or any available optional equipment without first completely understanding these instructions and any additional instructional material provided with this product. If you are unable to understand the warnings, cautions or instructions, contact a healthcare professional, dealer or technical personnel before attempting to use this equipment - otherwise injury or damage may occur.**
- **DO NOT** use substitute parts if components are damaged or missing. Contact your dealer immediately.
- Adjustable overbed tables have a powerful internal spring. Exercise extreme caution when working near the spring.
- **DO** regularly check frame screws to ensure they are tight.
- **DO** regularly check to ensure castors roll freely and locks work properly.
- The deluxe overbed table (model no. PB3037) has a weight limitation of 25 lbs.
- After any adjustments, repair or service and prior to use, be sure all attaching component parts are secure.

⚠ ACCESSORIES WARNING

Probasics products are specifically designed and manufactured for use in conjunction with Probasics accessories. Accessories designed by other manufacturers have not been tested by Probasics and are not recommended for use with Probasics products.

**For more information regarding Probasics products, please visit
www.pmiusa.biz**

Deluxe Overbed Table

Model No.PB3037

Overbed Table Assembly

WARNING:DO NOT remove the shipping screw until overbed table assembly is complete.

1. Position column upside down and insert release lever. (The release rod inside the column should be visible in the lower portion of the column. If not, shake the column).
 2. The curved side of the release lever should be closest to the top of the table and between the release rod and the column.
 3. Put table top upside down on a flat surface.
 4. Hold release lever in column and locate into top bracket---align holes---install bolts, first through top bracket, column and release lever, then second through top bracket and column. Tighten locknuts.
 5. Install a bolt through the bolt plate and column---hand tighten---repeat for other side.
 6. Put overbed table on its side to ease pressure on the bolts---align column and base and tighten both bolts.
 7. Remove the shipping screw from the column and replace with Phillips screw provided.
 8. The Overbed Table can now be adjusted to the desired height.
 9. To RAISE the table top, place your hand underneath the table top and move it upward to the desired height. Remove hand pressure to set height. Tighten the adjustment knob to secure the height of table.
 10. To LOWER the table top loosen the adjustment knob, lift the release lever and push the table top downward to the desired height. Let go of the release lever to set height. And re-tighten the adjustment knob.
- (The column locks in place when the release lever is released.)

