Knee Scooter

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OPERATING INSTRUCTIONS

Congratulations on the purchase of a Roscoe Knee Scooter. This medical device is designed to maximize mobility during your recovery from foot, ankle or leg surgery or injury.

Please make sure you read all instructions and warnings prior to attempting to assemble or operate this knee scooter. Otherwise, injury or damage may result.

PLATFORM REMOVAL

1. Remove knee scooter from box and remove all packaging including zip ties, bubble wrap and plastic on the wheels.
2. Carefully place knee scooter upside down so wheels are facing upwards.
3. Place two hands on knee platform and pull up to remove. If platform does not remove with slight effort, continue to pull and twist the knee platform left to right and repeat until platform becomes free.
4. Once the seat has been removed, carefully place the knee scooter on the ground so all four (4) wheels are flat.

STEERING COLUMN SETUP

Remove all packaging including zip ties, bubble wrap and plastic from the wheels before use.

The steering column will be in a folded position when unpacking it from the shipping carton. This feature also makes storing and transporting the scooter quick and convenient during everyday use.

1. To Raise the Steering Column

The quick release clamp ③ at the base of the steering column has an added security pin to prevent it from collapsing during use.

1. Release the lever and swing it down and away from the steering column.
2. Push the lever to the left to disengage the pin.
3. Lift the steering column until it's in the full upright position. Now, release the lever to re-engage the safety pin (the steering column will stand in the upright position if the safety pin is engaged).
4. Swing the lever back up and push it down so the end of the lever is pointing towards the ground to lock the steering column in the upright position.

To adjust the pressure on the lever, point the end of the lever up and rotate the lever counter clockwise to loosen and clockwise to tighten.

2. Handle Bar Height Adjustment

Handle bars should be set at waist level to keep your back straight. Stand on the side of the scooter where your pushing foot will be. Lock the parking brake (see step #2 of the braking section for details). Unscrew the height adjustment knob ② and adjust the height of the hand grips to a comfortable waist high fit. Insert and tighten the height adjustment knob (**DO NOT** over tighten).

STORAGE/TRANSPORTATION

The steering column can fold down for compact storage and transportation when not in use. Pull the lever up to release the clamp from the steering column and swing the lever forward to disengage. Push the lever to the left to disengage the pin and gently guide the steering column downward.

KNEE PLATFORM ADJUSTMENT

The knee platform can be used for right or left leg/foot injuries. The offset position of the platform provides better stability and increased foot clearance for pushing forward (always wear non-slip shoes on your non-injured foot). Look under the knee platform and follow the directional labels to know which pad should be facing forward for your type of injury.

1. Adjust the knee platform for a right or left leg injury by pulling out the post clamp lever ⑤ and removing the post pin ⑥. This will allow you to remove the platform to and turn it 180 degrees.
2. Reinsert the platform post and set the desired height. Once in place, insert the post pin and push down to lock the post clamp.
3. The vacant area between the pads is designed for the lip of a cast or boot. This area can be adjusted by loosening the four (4) screws under the knee platform and sliding to the desired position. This configuration is designed for maximum user comfort.

BRAKING

1. Applying The Braking

The braking wheel must be in contact with the floor, and the brake works best with weight on the knee scooter. If the brake fails to engage, the cable is probably too loose or too tight. Please see brake adjustment section.

- To use the brake, simply squeeze the lever with your finger tips toward the handle bar.
- To release the brake, simply release pressure on the brake lever.

2. Applying The Parking Brake

- To set the brake, push the brake lever in a downward direction away from the handlebar.
- To release the brake, squeeze the brake lever.

MOVEMENT

- Place your hands on the handlebars and apply the brake.
- Place the knee of your injured foot/leg on the front pad.
- Release the brake and propel yourself forward or backward with your non-injured foot, keeping the knee of your injured foot/leg on the pad (always wear non-slip shoes on your non-injured foot).
- Move slowly so that you can stop safely using your pushing foot.
- The hand brake may be used if needed but is primarily designed to maintain a stopped stance.

BRAKE ADJUSTMENT

The pressure of the brake assembly on the brake wheel should be enough so that when the hand lever is in the fully squeezed position, the brake wheel does not turn, but it is not overly tight. The tightness can be adjusted with the screws at either end of the brake cable.

- To correct, adjust cable with one or both screws at each end of the cable.
- Loosen the nut(s) and turn the adjustment screw(s) clockwise to loosen the cable or counterclockwise to tighten it. When adjustments are correct, the brake assembly will engage and the brake wheel will not turn when attempting to do so by hand.
- Return each nut to its position against the brake lever housing and cable mounting bracket at the lower end of the cable.

If these adjustments do not correct the problem, it may be necessary to make additional adjustments where the brake cable is attached to the brake shoe.

- Loosen the screw and nut that holds the brake cable; located on the rear brake assembly ⑦.
- Inspect the cable for fraying and replace if needed.
- Set adjustment screws at each end of the cable in the mid-position.
- Pull the brake cable with a pair of pliers to remove any slack in the cable and re-tighten the screw and nut.

NOTE: Bicycle mechanics are a good source for advice, repairs and service.

WARNINGS

DO NOT ATTEMPT TO SET-UP, ADJUST OR OPERATE THIS KNEE SCOOTER WITHOUT READING ALL INSTRUCTIONS CAREFULLY.

- As with any new physical activity, instruction, practice and common sense are necessary for successful use of this scooter. Your health professional or Roscoe Knee Scooter distributor can provide instruction and a demonstration. Practice operating this scooter to increase your skill and confidence. Gradually increase the time you spend using the scooter to allow your muscles to adjust to this new activity.
- Failure to properly engage the quick-release clamp may cause the steering column to move while the scooter is in use and may cause you to lose control. With the quick release clamp properly adjusted, the handle bars will be locked in place.
- Be sure to check the quick-release clamp is secure and tight before each use or injury may result.
- Injury may result if seat knee platform adjustment instructions are not followed.
- **DO NOT** use the Roscoe Knee Scooter to pull yourself up
- The scooter is intended for use with a single person in place on the scooter. ALL WHEELS must be in contact with the floor/ground at ALL TIMES. Failure to maintain this could result in improper balance of the scooter and cause injury or damage.
- **ALWAYS INSPECT SCOOTER PRIOR TO RIDING.** Properly inspecting and maintaining your scooter can reduce the risk of injury. Make sure that all parts are secure and are in good working order. Always check the fasteners, wheels, brakes and seat for tightness prior to use. If scooter is not in good working order or fasteners cannot be tightened, **DISCONTINUE USE OF THE SCOOTER.**
- This is not a toy. **DO NOT** let children play on scooter.
- Always wear protective footwear, never ride barefooted or in sandal, and keep shoelaces, loose clothing and/or accessories out of the way of the wheels and handle bars.

- Be aware of your surroundings when operating the scooter. Look for hazards and avoid them. Avoid streets and surfaces with water, sand, gravel, dirt, leaves, and other debris. Wet, slick, uneven and/or rough surfaces may impair traction and contribute to possible accidents. **DO NOT** ride a scooter in mud, ice, puddles or water.
- Never use near steps, sloped driveways/hills with an incline of more than 15 degrees or swimming pool areas.
- Move slowly so that you can stop safely using your pushing foot.
- **DO NOT** use the hand brake as a primary stopping tool. It may be used if needed but is primarily designed to maintain a stopped stance.
- **DO NOT** use while under the influence of mind-altering drugs, alcohol or while suffering from dizziness.
- Use caution when moving from carpeted to hard surface floors or cracks/separations in walkway surfaces; avoid sharp bumps, drainage grates, and sudden surface changes. Scooter may abruptly stop (i.e. moving from a curb to the street and back onto the sidewalk).
- Never make a sharp turn at high speed; do not turn at speeds over 1 mph (very slow walking speed). **DO NOT** turn unless one foot is on the ground.
- **DO NOT** sit on or rest knee on the scooter while it is being loaded into or being transported in a vehicle.
- Always keep both hands on the handle grip when scooter is in use.
- **DO NOT** attempt to reach objects that are out of your immediate reach while using this device. **DO NOT** reach for objects if you have to lean or shift position on the scooter. **DO NOT** lean forward, backward or to either side while using the scooter.
- **DO NOT** store in direct sunlight or outdoors. Exposure to extreme temperatures (heat, cold, moisture, etc.) may make surfaces (knee platform, hand grips, etc.) hot to the touch, weak or brittle. Ensure hand grips **DO NOT** twist and skin contact surfaces are cool to the touch before use.

FAILURE TO USE COMMON SENSE AND HEED THE ABOVE WARNINGS FURTHER INCREASES RISK OF SERIOUS INJURY. USE AT YOUR OWN RISK AND WITH APPROPRIATE AND SERIOUS ATTENTION TO SAFE OPERATION. USE CAUTION.

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WARRANTY

This warranty is solely for the benefit of the original consumer purchaser.

This device was built to exacting standards and carefully inspected prior to shipment. Your Roscoe brand product carries a limited lifetime warranty on the frame and 90 days on all other non-durable components including but not limited to rubber accessories, wheels, grips, brake cables, knee pads and seat, which are subject to normal wear and need periodic replacement. If within such warranty period, the product proves defective, the product will be repaired or replaced, at the option of Roscoe Medical. This warranty does not cover device failure due to owner misuse or negligence, or normal wear and tear.

This warranty shall not apply to serial numbered products if the serial number has been removed or defaced, or to defects or damages caused by unauthorized parts, attachments or repairs, negligence, accident, or any circumstance beyond Roscoe Medical's control, such as improper operation, maintenance or storage.

If you have a question about your Roscoe product or this warranty, please contact an authorized Roscoe Dealer.

Roscoe Medical
6753 Engle Road
Middleburg Heights, Ohio 44130
Ph: (800) 3-ROSCOE (376-7263)
www.roscoemedical.com
42-ROS-KSB_03 ©2016