

Patient Interfaces

Roscoe Full Face Mask

Roscoe Medical's Full Face Mask with Comfort Cushion™ is made of high-quality clear silicone with a virtually indestructible frame. The lightweight mask comes complete with headgear and dual 360 degree swivel at the tubing attachment point. Mask is dishwasher safe and comes complete with headgear.



Product Features

- Made of high-quality, lightweight, clear silicone with a flexible frame
- Dual 360° swivel at the tubing attachment point
- Made with half the parts of popular brands, meaning less patient confusion
- Headgear with Velcro® fasteners provide maximum versatility and movement
- Latex-free
- Dishwasher safe for easy cleaning
- Full face mask comes complete with headgear (black)
- 90 day limited warranty
- HCPCS Code - A7030* (mask)
- HCPCS Code - A7035 (replacement headgear)



Roscoe Medical Offers a 30 Day Money-Back Guarantee on all Roscoe® CPAP Masks, effective from the date of purchase.

If your patient is not 100% completely satisfied and it is within 30 days from your Roscoe purchase date, we will refund the purchase price to you, no questions asked.

Item #	Description	Qty
CPM-FFS	Full Face Mask with Headgear, small	each
CPM-FFM	Full Face Mask with Headgear, medium	each
CPM-FFL	Full Face Mask with Headgear, large	each
CPH-FF	Replacement Headgear, black	each

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



To Order:
 (800) 3-ROSCOE
 fax: (440) 572-4261
 roscoemedical.com