

# Electrotherapy

## InTENSity™ 5000 Hybrid

The InTENSity™ 5000 Hybrid device combines both digital and analog technology to create an innovative pain management device that is easy-to-use. The user will receive the right therapy at the right time by using traditional analog dials to manage intensity levels through LED technology. The InTENSity™ 5000 Hybrid comes complete with five modes, a timer, electrodes, lead wires, carrying case, user manual and 9-volt battery.



### Product Features:

- Both digital and analog unit
- Comes complete with carrying case, two sets of lead wires, 9-volt battery, pack of four self-adhesive reusable electrodes and manual
- Unit includes five modes of operation - Burst, Normal, Modulation 1, Modulation 2 and SD mode
- Dimensions: 4.6" x 2.6" x 11.2"
- Weight: 0.2 lb. without battery
- One year warranty
- HCPCS Code: E0730\* for TENS Unit
- HCPCS Code: A4595\* for electrotherapy supplies



### Product Specifications:

- Power Source: 9-volt battery or AC adaptor (optional)
- Channel: Dual isolated between channels
- Pulse Amplitude: Adjustable 0-80 mA, peak into 500 ohm load each channel
- Pulse Rate: 2-150 Hz
- Burst Rate: 0-5 Hz (Burst mode)
- Pulse Width: 30-260 microseconds
- Wave Form: monophasic square pulse wave
- Timer: 15, 30 and 60 minutes or continuous

Item #	Description	Qty
DI3302	InTENSity™ 5000 Hybrid	each
DI0001X	InTENSity™ Series AC Adaptor	each

\* Prescription Required: This is a Class II Medical Device. Federal law restricts these devices to sale by or on the order of a licensed health care professional.

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



**To Order:**

(800) 3-ROSCOE  
 fax: (440) 572-4261  
 roscoemedical.com