

SOUNDCARE™ PLUS

PROFESSIONAL ULTRASOUND DEVICE

Unique in nature, the SoundCare™ Plus professional ultrasound device comes complete with two sound heads that can be plugged in simultaneously for more targeted therapy. Simple to operate, the user can switch between wands with the single push of a button.

FEATURES:

- Comes complete with two sound heads, 1 cm and 5 cm for more targeted therapy
- Used for treatment of chronic and acute muscular pain
- User can switch between dual soundheads with the push of a single button
- Fully functional 1 MHz and 3 MHz frequencies
- 20 user defined presets with pulsed and continuous therapy
- Ergonomic handle design
- Two year warranty

ORDER DETAILS:

DQ9275	SoundCare™ Plus Professional Ultrasound Devices
LC2828	TheraSonic™ Ultrasound Gel, 2 ounces
LS5288	TheraSonic™ Ultrasound Gel, 8 ounces
LS5266	TheraSonic™ Ultrasound Gel, 16 ounces

*Prescription Required: This is a Class II Medical Device. Federal law restricts these devices to sale by or on the order of a licensed health care professional.



TECHNICAL SPECIFICATIONS

- **Acoustic Frequency:** 1 MHz \pm 10%, 3 MHz \pm 10%
- **Output Power:** 0.5 - 10.0 W \pm 20%, when duty factor \geq 80% for 5 cm²
0.5 - 15.0 W \pm 20%, when duty factor \leq 70% for 5 cm²
0.1 - 2.0 W \pm 20%, when duty factor \geq 80% for 1 cm²
0.1 - 3.0 W \pm 20%, when duty factor \leq 70% for 1 cm²
- **Pulse Repetition Rate:** 100 Hz \pm 10%
- **Duty Factor:** 10% - 100%, Stepping 10%
- **Timer:** Max 30 minutes
- **Effective Radiating Area:** 5.0 cm² \pm 20%; 1.0 cm² \pm 20%
- **Actual Intensity:** 3.0 w/cm² \pm 20% (1.0 MHz)
3.0 w/cm² \pm 20% (3.0 MHz)
- **Raw (Max):** 5.0
- **Beam Type:** Collimated
- **Material of Applicator:** Aluminum
- **Waterproof Grade:** IPX7 only for treatment head