

THERATOUGH™ UX2

Advanced Ultrasound Device

Equipped with a menu-driven, color touchscreen interface that intuitively groups and displays clinical protocols, guiding clinicians step-by-step to ensure the ideal therapeutic treatment.

- Appropriate parameter settings
- Ultrasound applicator recommendations
- Dual-holster design to accommodate both size applicators
- Ergonomic designed sound head provides for an excellent ERA

A few simple key presses are all that is needed to quick-start a treatment. Parameters can be easily selected and adjusted. Clinicians can choose to use preset clinical protocols or customize their own. The “favorites” feature allows user-defined protocols to be stored for future use.

The UX2 offers a comprehensive library with treatment education for better treatment outcomes and improved clinic efficiency.

FEATURES:

- Comprehensive Library Provides Treatment Education
- Sleek, User-Friendly, Ergonomic Design for Tabletop or Therapy Cart Use
- Large, Easy-to-Read, Color LCD Display with Touch Screen Interface
- Clinical Protocols Provide Parameters for a Variety of Acute and Chronic Conditions
- “Favorites” Feature Allows Storage Up to 99 Custom Programs for Future Use
- 1 and 3 MHz Ultrasound Frequencies
- Pulsed and Continuous Modes
- Comes with 1 cm² and 5 cm² Applicators
- Low BNR Ratio 5:1
- Ultrasound Amplitude Can Be Toggled Between W and W/cm²
- Sound Head Coupling Feature
- Optional Therapy Cart Integrates with the Unit and Keeps Supplies Organized and Ready for Treatments



Optional
Therapy
Cart



TECHNICAL SPECIFICATIONS

ULTRASOUND

- **Frequency:** 1 and 3 MHz
- **Duty Cycles:** Pulsed and Continuous
- **Pulse Frequency:** 100 Hz
- **Pulse Duration:** 1-9 mSec ± 10% (set by Duty Cycle)
- **Ramp (Fixed):** 1s
- **Treatment Time:** 1 - 60 minutes

ORDER DETAILS:

DQ2000	TheraTouch™ UX2
DQCART-2	Therapy Cart with Storage