



Bariatric Commode Extra Wide Seat

ProBasics Bariatric Commode Extra Wide Seat features a 24-inch wide platform seat and reinforced steel frame to support patient weights up to 650 pounds and is ideal for larger patients. It comes complete with commode bucket, splash guard and lid. The versatile design includes a commode bucket with cover, armrests for added comfort, angled legs for stability and non-marring rubber tips.

FEATURES

- Welded steel construction with 24-inch extra wide snap on platform seat designed for strength and stability
- Tool-free assembly
- Commode opening is contoured for increased comfort
- Seat height readily adjusts in one-inch increments
- Angled legs for stability with non-marring rubber tips
- Plastic armrests provide added comfort
- Comes complete with a commode bucket and cover
- Color: gray frame, gray accessories
- HCPCS Code: E0168

SPECIFICATIONS

Width b/w Arms:	26"
Outside Width:	27"
Seat Dimensions:	24"W x 18"D
Seat Opening:	9"W x 10.5"D
Adjustable Seat Height:	17" – 23"
Color:	Gray
Weight Capacity:	650 lb
Carton Dimensions:	32.5" x 22" x 26"
Shipping Weight:	41 lb
Warranty:	Limited lifetime on frame One year on parts

Item #	Description	UPC	UOM
BSB24C	ProBasics Bariatric Commode with Extra Wide Seat	815067071890	2/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.