


4580P



4580P	DESCRIPTION	ITEM #
1	Mattress Only	3102M
2	Top Cover	C_3102C
3	Base Cover	PRT-3102-BASE
4	Inside Base Cover w/Cells	PRT-3102-IBCT
5	LAL Cell	PRT-3102-VCCELL
6	Non-LAL Cell	PRT-3102-NVCELL
7	CPR Plug	PRT-CPR-RIBBON
8	Base Cover (not shown)	C_3102-B-88
9	Pump	4253E
10	Mattress Hose	90165
11	Check valve (not shown)	PRT-CHK-VALVE

Images are for reference purposes only and may not resemble the specific APM model.