

ARUBA-42

BARIATRIC GROUP 1 FOAM MATTRESS 600 LB WEIGHT CAPACITY, 42" X 80" X 6"

*Ships standard
with select bariatric
bed packages!*


The Meridian Medical Aruba-42 joins the Compass Health line of therapeutic Group 1 support surfaces. This bariatric layered foam mattress arrives compressed for easy handling and delivery, plus the smaller carton (just 14" x 14" x 44 ½") means lower shipping costs!

Multi-density layered foam provides a surface that helps reduce the occurrence of pressure sores and the heel foam section aids in preventing pressure sores on patients' feet. With a urethane-backed, water-resistant zippered cover, cleanup is quick and easy.

In just a few minutes the mattress is ready for patients up to 600 lbs. to enjoy a comfortable 42" x 80" x 6" sleeping surface.

FEATURES:


- Boxed and compressed for easy handling & shipping savings (just 14" x 14" x 44 ½")
- Group 1 multilayer foam helps reduce and prevent pressure sores
- Supports patients up to 600 lbs.
- Fully decompresses in just a few minutes to 42" x 80" x 6"
- Water-resistant zippered cover for easier cleaning
- Includes a protective fire barrier (fire sock)
- 1-year limited warranty



Foam colors may vary

PRODUCT SPECIFICATIONS:

Shipping Dimensions:	14" x 14" x 44 ½", 27.6 lbs.
Dimensions:	42" x 80" x 6"
Weight Capacity:	600 lbs
Product Weight:	24.25 lbs
Foam Density:	1.6 PCF / 1.85 PCF
Warranty:	12-Month Limited
Safety:	Meets 16 CFR Part 1633, No flame retardant chemicals
Prop 65:	Yes (DEHP)



Give patients the necessary additional comfort and support to prevent pressure sores, choose the ARUBA-42 Foam Mattress.

Item #	Description	Weight Capacity	Qty
ARUBA-42	42" Group 1 Foam Mattress	600 lbs	each