

## BreatheRite™ Rigid Valved Holding Chamber and Masks



The BreatheRite™ Valved Holding Chamber (VHC) combines design advantages with portability to provide both effective medication delivery and encourages increased patient use.

### PRODUCT SPECIFICATIONS

- Anti-static chamber enables effective medication delivery
- Designed to promote greater drug volume output
- Improper inspiratory rate is indicated by a soft vibrating sound
- Chamber design eliminates the need for a flow whistle
- Face mask (all detachable) options available; can be used with or without mask



Item #	Description	Case Quantity	UOM
BR300	BreatheRite Aerosol Chamber	20	case
M10	Adult Mask	10	case
M20	Pediatric Mask	10	case