



## Three-in-One Steel Commode

The **ProBasics Three-in-One Steel Commode** has three versatile functions; it can be used as a toilet safety frame, a commode, and a raised toilet seat. It features a removable back and comes complete with commode bucket, splash guard and lid. The versatile design includes an elongated deeper seat, armrests for added comfort, angled legs for stability and non-marring rubber tips.

### FEATURES

- Can function as a free-standing commode, a toilet safety frame (with backrest removed), or as a raised toilet seat with arms
- Features deeper elongated seat depth for less mess and to allow more room for cleaning
- Durable steel-welded frame is attractive, strong and easy-to-clean
- Tool-free assembly, includes plastic snap on seat with lid
- Seat height readily adjusts in one-inch increments
- Plastic arms provide added comfort
- Comes complete with commode bucket with handle, cover and splash guard
- HCPCS Code: E0163\*

### SPECIFICATIONS

Width between arms:	17.5"
Outside Width:	24"
Seat Dimensions:	13.5"W x 16.5"D
Seat Opening:	8"W x 10"D
Adjustable Seat Height:	15.5" – 23.5"
Color:	Gray
Weight Capacity:	350 lb
Carton Dimensions:	25.5" x 25.5" x 25"
Shipping Weight:	57 lb
Warranty:	Limited lifetime on frame One year on parts

Item #	Description	UPC	UOM
BS31C	ProBasics Three-in-One Commode with Elongated Seat	815067071838	4/cs

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.