



Transfer Bench

The **ProBasics Transfer Bench** helps make transfers in and out of the bathtub safer and more comfortable. Constructed with an aluminum anodized frame and blow-molded bench, it features a reversible set back for right or left handed transfers. The angled legs have one-inch, easy-snap button adjustments and wide-based, non-skid feet for added sturdiness and safety. Tool-free assembly and limited lifetime warranty are standard.

FEATURES

- Aluminum frame is lightweight, sturdy and corrosion-proof
- Durable blow-molded bench and back have a textured surface and are contoured with drainage holes to reduce slipping
- Wide-based tips and suction cups add sturdiness and prevent slipping
- Legs have one-inch easy snap-button adjustments
- Converts easily for left or right hand transfers
- Quick tool-free assembly
- HCPCS Code: E0247*

SPECIFICATIONS

Overall Footprint:	28.5"W x 17.5"D
Seat Dimensions:	25"W x 16"D
Adjustable Seat Height:	19" - 23"
Weight Capacity:	300 lb
Carton Dimensions:	28" x 17" x 11.5"
Shipping Weight:	23 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	UPC	UOM
BSTB	ProBasics Transfer Bench	815067071968	2/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.