

Respiratory – Aerosol Therapy



BEAR PEDIATRIC NEBULIZER SYSTEM

Designed to reduce anxiety during treatments and encourage compliance, the child-friendly Bear Compressor Nebulizer provides effective aerosol treatments with a state-of-the-art piston compressor. Delivers 30 PSI at maximum pressure and 10 PSI operating pressure (liter flow range 5 - 8 LPM). Equipped with a convenient nebulizer parking post. Comes complete with compressor, disposable nebulizer, angled mouthpiece, 7-foot kink-resistant tubing, filter, and a five year warranty.



Specifications:

Method of Nebulization:	Piston Compressor
Medication Capacity:	2 - 8 ml
Particle Size:	0.5 µm to 10 µm
Unit Dimensions:	7.15" (L) x 6.90" (W) x 5.5" (H)
Product Weight:	3.4 lbs.
Maximum Pressure:	30 PSI
Operating Pressure:	10 PSI
Liter Flow:	8 LPM
Electrical Requirement:	120 VAC/60Hz
Warranty:	Limited 5 Year

Features:

- Child-friendly design helps to reduce anxiety during treatments
- Provides smooth, quiet treatments
- Equipped with a convenient nebulizer parking post
- Complete with compressor, nebulizer, angled mouthpiece, tubing, filter and instruction manual
- Nylon carry bag available: BAG-DRW

HPCPS Codes:

E0570*	Nebulizer with Compressor
A7003*	Disposable Nebulizer Kits
A7005*	Reusable Nebulizer Kits
A7013*	Foam Filters
A7015*	Aerosol Mask

Item #	Description	UOM
JLN-2315AS	Bear Pediatric Nebulizer	each
Accessory		
BAG-DRW	Nylon Carry Bag for Pediatric Neb	each

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

** Prescription Required: This is a Class II Medical Device. Federal law restricts these devices to sale by or on the order of a licensed health care professional.