

MICiTECH



10 Liter Stationary Oxygen Concentrator

Introducing the MICiTECH 10 Liter Stationary Oxygen Concentrator. This quiet, compact and durable 10 Liter concentrator is low maintenance and an affordable choice for Home Medical Equipment and Hospice Dealers.

PRODUCT FEATURES:

- Compact and easy to transport
- Quiet operation
- Durable double-wobble compressor
- Designed for 24/7 use
- Low maintenance
- Easy removal and replacement of cabinet and compressor intake filters
- 3-Year Limited Warranty

TECHNICAL SPECIFICATIONS:

Product Name:	10 Liter Oxygen Concentrator	Electrical Requirements:	AC110V / 60 Hz
Dimensions:	15.74"L x 11.81"W x 21.26"H	Power Consumption:	700W
Weight:	44.6 lbs	Rated Current Input:	7.0A
Liter Flow:	0.5 -5 LPM	Sound Level:	≤55 dBA
Flow Rates:	1.0 ~ 10 LPM, increments of 1 LPM	Operating Environment:	Temperature: 41° ~ 104°F Humidity: </= 80% Atmospheric Pressure: 12.47 ~ 15.37 PSI
Maximum Outlet Pressure:	8.7 ~ 11.6 psi	Shipping/Storage:	Temperature: -4° ~ 131°F Humidity: 93% Atmospheric Pressure: 7.25 ~ 15.37 PSI
Oxygen Concentration:	93% (±3%)	Standards:	AAMI ANSI ES60601-1 IEC 60601-1-2 IEC60601-1-8 IEC 60601-1-11 ISO 80601-2-69
Mode of Operation:	Continuous Flow	Warranty:	3 Year Limited

Item #	Description	UPC	Qty
CP101	MICiTECH 10 Liter Stationary Oxygen Concentrator	860012953332	each
X1MMCP523	Compressor Intake Filter	N/A	each
X1MMCP505	Cabinet Filter	N/A	20/pk
ILB-516S	In-Line Bacteria Filter	N/A	each