

MICiTECH



5 Liter Stationary Oxygen Concentrator

Introducing the MICiTECH 5 Liter Stationary Oxygen Concentrator. This quiet, compact and durable 5 Liter concentrator is low maintenance and an affordable choice for Home Medical Equipment and Hospice Dealers.

PRODUCT FEATURES:

- Compact and easy to transport
- Quiet operation
- Durable double-wobble compressor
- Designed for 24/7 use
- Low maintenance
- Easy removal and replacement of cabinet and compressor intake filters
- 3-Year Limited Warranty

TECHNICAL SPECIFICATIONS:

Product Name:	5 Liter Oxygen Concentrator	Electrical Requirements:	AC110V / 60 Hz
Dimensions:	16.04"L x 11.81"W x 21.26"H	Power Consumption:	390W
Weight:	32 lbs	Rated Current Input:	3.2A
Liter Flow:	0.5 -5 LPM	Sound Level:	≤50 dBA
Flow Rates:	0.5 ~ 5 LPM, increments of 0.5 LPM	Operating Environment:	Temperature: 41° ~ 104°F Humidity: </= 80% Atmospheric Pressure: 12.47 ~ 15.37 PSI
Maximum Outlet Pressure:	5.8 ~ 8.7 psi	Shipping/Storage:	Temperature: -4° ~ 131°F Humidity: 93% Atmospheric Pressure: 7.25 ~ 15.37 PSI
Oxygen Concentration:	93% (±3%)	Standards:	AAMI ANSI ES60601-1 IEC 60601-1-2 IEC60601-1-8 IEC 60601-1-11 ISO 80601-2-69
Mode of Operation:	Continuous Flow	Warranty:	3 Year Limited

Item #	Description	UPC	Qty
CP502	MICiTECH Stationary 5 Liter Oxygen Concentrator	860012953318	each
X1MMCP504	Compressor Intake Filter	N/A	each
X1MMCP505	Cabinet Filter	N/A	20/pk
ILB-516S	In-Line Bacteria Filter	N/A	each