



Medical Gas Cylinders

Compass Health offers **empty aluminum oxygen cylinders** in standard sizes packaged to meet your requirements. You can order units in packs of six or in competitively priced pallet quantities.



NOTE: CYLINDER SIZE SHOWN IS NOT PROPORTIONALLY CORRECT

- All cylinders are available in both CGA-870 wrench and toggle valves (excluding the CGA-540 configured sizes)
- Dual control valve option available upon request
- Cylinder finish is easy to clean and maintain
- Private labeling options available with your choice of stamping or silk screen. Pallet minimums are required, additional fees may apply.

Item #	Cylinder Size	Dimensions	Oxygen Capacity (Liters)	Valve Type	Qty/UOM
PX-8704-1T	M4	9" x 3"	113 Liters	Toggle	6/pack
PX-8704-1TB	M4	9" x 3"	113 Liters	Toggle	144/pallet
PX-8704-1W	M4	9" x 3"	113 Liters	Wrench	6/pack
PX-8704-1WB	M4	9" x 3"	113 Liters	Wrench	144/pallet
PX-8700-1T	M6	11 3/4" x 3"	164 Liters	Toggle	6/pack
PX-8700-1TB	M6	11 3/4" x 3"	164 Liters	Toggle	144/pallet
PX-8700-1W	M6	11 3/4" x 3"	164 Liters	Wrench	6/pack
PX-8700-1WB	M6	11 3/4" x 3"	164 Liters	Wrench	144/pallet
PX-8705-1T	ML6	7 3/4" x 4 1/2"	170 Liters	Toggle	6/pack
PX-8705-1TB	ML6	7 3/4" x 4 1/2"	170 Liters	Toggle	81/pallet
PX-8705-1W	ML6	7 3/4" x 4 1/2"	170 Liters	Wrench	6/pack
PX-8705-1WB	ML6	7 3/4" x 4 1/2"	170 Liters	Wrench	81/pallet
PX-8701-1T	M9	10 1/2" x 4 1/2"	246 Liters	Toggle	6/pack
PX-8701-1TB	M9	10 1/2" x 4 1/2"	246 Liters	Toggle	81/pallet
PX-8701-1W	M9	10 1/2" x 4 1/2"	246 Liters	Wrench	6/pack
PX-8701-1WB	M9	10 1/2" x 4 1/2"	246 Liters	Wrench	81/pallet
PX-8702-1T	D	16" x 4 1/2"	425 Liters	Toggle	6/pack
PX-8702-1TB	D	16" x 4 1/2"	425 Liters	Toggle	81/pallet
PX-8702-1W	D	16" x 4 1/2"	425 Liters	Wrench	6/pack
PX-8702-1WB	D	16" x 4 1/2"	425 Liters	Wrench	81/pallet
PX-8703-1T	E	25" x 4 1/2"	679 Liters	Toggle	6/pack
PX-8703-1TB	E	25" x 4 1/2"	679 Liters	Toggle	81/pallet
PX-8703-1W	E	25" x 4 1/2"	679 Liters	Wrench	6/pack
PX-8703-1WB	E	25" x 4 1/2"	679 Liters	Wrench	81/pallet



Medical Gas Cylinders

HOMEFILL-STYLE OXYGEN CYLINDERS

Compass Health offers empty aluminum oxygen cylinders with **Invacare® HomeFill®-Style Valves** in four sizes packaged to meet your requirements. Available in single pack or 6 pack quantities.

- CGA-870 wrench valve available
- Brushed aluminum with clear coat finish
- Ships empty
- For use with Invacare HomeFill units
- Free freight on 6 packs for orders shipped within the Continental United States only

Item #	Description	UOM	Qty
353763	Invacare HomeFill-Style ML6 Cylinder with Post Valve and Transfill Port	Pack	6
HFML6-EA	Invacare HomeFill-Style ML6 Cylinder with Post Valve and Transfill Port	Each	1
353772	Invacare HomeFill-Style M9/C Cylinder with Post Valve and Transfill Port	Pack	6
HFC-EA	Invacare HomeFill-Style M9/C Cylinder with Post Valve and Transfill Port	Each	1
353750	Invacare HomeFill-Style D Cylinder with Post Valve and Transfill Port	Pack	6
HFD-EA	Invacare HomeFill-Style D Cylinder with Post Valve and Transfill Port	Each	1
353769	Invacare HomeFill-Style E Cylinder with Post Valve and Transfill Port	Pack	6
HFE-EA	Invacare HomeFill-Style E Cylinder with Post Valve and Transfill Port	Each	1

CYLINDER ACCESSORIES



PV-TG



PV-WR



ET-CW



2544-2



90003



SMK-SGN

Item #	Description	Qty
90003	Replacement T handle for all Roscoe Regulators	each
PV-TG	Post valve, toggle-style, CGA-870	each
PV-WR	Post valve, wrench style, CGA-870	each
ET-CW	Cylinder wrench for wrench valve	25/pk
MPV-4046	Toggle replacement, metal	each
MPV-4047	Replacement Pin	each
MCW-2-1	Cylinder wrench chain with clip	each
2544-2	Regulator yoke washer with rubber ring	25/pk
47012	Medical post valve locking band, white	each
SMK-SGN	"No Smoking, Oxygen in Use" sign, non-adhesive, printed in red	100/pk

Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Compass Health Brands unless otherwise noted.