



INTENSITY™ 10

Portable Electrotherapy

The InTENSity 10 Digital TENS machine allows patients to select from 10 preset treatment protocols for targeted TENS electrotherapy. With the press of two buttons, this easy-to-use TENS device is ready to operate effectively based on the area of the body that is being treated. Comes complete with AC adapter and 9-volt battery to deliver effective therapy on the go.

FEATURES:

- Simple to Operate, Portable TENS Unit for Pain Management Can Be Used in the Comfort of Your Own Home
- Device Has Ten Preset Body Part Indications to Select From, Making It Easy for Patients to Target Specific Areas for Pain Relief
- Body Part Selections Include: Neck, Shoulder, Mid Back, Low Back, Elbow, Wrist, Hip, Knee, Ankle, and Foot
- Each Individual Body Part Selection Has a Customized Treatment Program Based on Researched Clinical Treatment Protocols
- Portable TENS Machine Comes Complete with Carrying Case, Two Sets of Lead Wires, 9-Volt Battery, AC Adapter, Manual, and Pack of Four Self-Adhesive Reusable Electrodes
- HCPCS Codes: E0730, A4595
- One-Year Warranty
- **Prescription Required:** This Is a Class II Medical Device. Federal Law Restricts These Devices to Sale By or on the Order of a Licensed Health Care Professional



TECHNICAL SPECIFICATIONS

- **Power Source:** 9V Battery, AC Adapter
- **Dimensions:** 4.5" x 2.55" x 0.9"
- **Weight:** 0.28 lb
- **Pulse Amplitude:** Adjustable 0-105mA at 1000 ohm load each channel
- **Pulse Width:** 100 to 260µs
- **Pulse Rate:** 50-150Hz
- **Waveform:** Bi-phase square pulse wave
- **Timer:** 15, 30, 60 minutes and continuous

ORDER DETAILS:

DI1010

InTENSity 10

