



INTENSITY™ 5000

Hybrid TENS

The InTENSity 5000 Hybrid combines both digital and analog technology to create an innovative pain management device that is easy-to-use. The user can adjust the unit by using the traditional analog dials to manage intensity levels through an LED display common with digital units. The InTENSity 5000 is equipped with five treatment modes and a timer.



FEATURES:

- Comes Complete with Carrying Case, Two Sets of Lead Wires, 9-Volt Battery, Pack of Four Self-Adhesive Reusable Electrodes and Manual
- Five Modes of Operation: Burst, Normal, Modulation 1, Modulation 2 and SD Mode
- HCPCS Code: E0730, A4595
- One Year Warranty
- Prescription Required: This Is a Class II Medical Device. Federal Law Restricts These Devices to Sale By or on the Order of a Licensed Health Care Professional.

ORDER DETAILS:

DI3302

Intensity 5000 – Hybrid TENS

TECHNICAL SPECIFICATIONS

- **Power Source:** 9-volt battery or AC adaptor (optional)
- **Channel:** Dual isolated between channels
- **Pulse Amplitude:** Adjustable 0-80 mA, peak into 500 ohm load each channel
- **Pulse Rate:** 2-150 Hz
- **Burst Rate:** 0-5 Hz (Burst mode)
- **Pulse Width:** 30-260 microseconds
- **Wave Form:** monophasic square pulse wave
- **Timer:** 15, 30 and 60 minutes or continuous
- **Dimensions:** 4.6" x 2.6" x 11.2"
- **Weight:** 0.2 lb. without battery