

INTENSITY™ TWIN STIM III

Portable Electrotherapy

The InTENSity Twin Stim III is a portable, digital device that combines TENS for pain management and EMS for muscle rehabilitation and re-education. The device comes equipped with four modes for TENS therapy and three modes for EMS therapy. The unit is among the upper echelon of pain relief devices, with a max output of 105mA. It also comes with an AC adapter and removable belt clip for easy at home and on the go use.



FEATURES:

- Intended for Use of Pain Relief, Muscle Reeducation, Increasing Local Blood Circulation and Maintaining Range of Motion
- Unit Combines TENS and EMS Therapy for Optimal Treatment
- Simple-to-Operate, Portable Electrotherapy Device for Pain Management
- Four Modes for TENS Therapy: Burst, Normal, Mod 1, and Mod 2
- Three Modes for EMS Therapy: Synchronous, Alternate, and Delay
- Comes Complete with Carrying Case, Pair of Lead Wires, 9-Volt Battery, AC Adapter, Manual, and Pack of Four Self-Adhesive Reusable Electrodes
- One-Year Warranty
- Rx Only
- HCPCS Code: E0730

TECHNICAL SPECIFICATIONS

- **Power:** 9V Alkaline DC
- **Waveforms:** TENS, EMS

	TENS	EMS
Wave Form	Monophasic Square	Biphasic
Pulse Rate	1 - 150 Hz	1 - 150 Hz
Pulse Width	50-300 μ s	50-300 μ s
Output Voltage	Constant Current	Constant Current
Current	0-105 mA	0-105 mA
Modes	4	3
Timer	1-60 Constant	1-60 Constant

ORDER DETAILS:

DI3717 InTENSity Twin Stim III