



INTENSITY™ IF COMBO

Portable Electrotherapy

The InTENSity IF Combo is a portable TENS unit that combines TENS and Interferential waveforms, giving you the best combination of treatment options for patients needing options beyond TENS alone. The IF Combo has four modes of operation in both TENS and Interferential, ensuring your patient is getting the most out their treatment options. The IF Combo has an output of up to 105 mA from a standard 9-volt battery, or with an AC adapter, so the patient can get relief virtually anywhere.



FEATURES:

- Simple-to-Operate, Portable Electrotherapy Device for Pain Management
- Combines TENS (Transcutaneous Electrical Nerve Stimulation) and Interferential (IFC) Treatment in a Single TENS Machine
- Equipped with Four Treatment Modes
- Comes Complete with Carrying Case, Pair of Lead Wires, 9-Volt Battery, AC Adapter, Manual, and Pack of Four Self-Adhesive Reusable Electrodes
- HCPCS Code: E0730 for TENS Unit
- HCPCS Code: E4595 for Electrotherapy Supplies
- One-Year Warranty
- Rx Only

TECHNICAL SPECIFICATIONS

- **Power Source:** Standard 9-volt battery or AC adapter
- **Dimensions:** 4.5" x 2.55" x 0.9"
- **Weight:** 0.28 lb

	TENS	Interferential
Wave Form	Monophasic Square	Biphasic Square
Pulse Rate	1 - 150 Hz	4kHz
Pulse Width	50-300 μ s	125 μ s: phase duration; 250 μ s
Output Voltage	Constant Current	Constant Current
Current	0-105 mA	0-105 mA
Modes	4	4
Timer	1-60 Constant	1-60 Constant, Freq: 1-50 beats

ORDER DETAILS:

DI4738

InTENSity™ IF Combo

