



INTENSITY™ SELECT COMBO

Portable Electrotherapy

The InTENSity Select Combo is a powerful, easy to use portable electrotherapy device. This unit features four therapeutic modes including TENS, EMS, IF and Microcurrent. These four modes and a strong current of 105mA provide a wide range of treatment options and flexibility to treat a variety of musculoskeletal issues. The InTENSity Select Combo can help patients with pain management, muscle re-education, muscle range of motion and more.

FEATURES:

- Four Waveforms: TENS (Transcutaneous Electrical Nerve Stimulation), EMS (Electrical Muscle Stimulation), IFC (Interferential) and Microcurrent
- Over 14 Programs to Choose From
- InTENSity Select Combo Tens Machine Comes with Carrying Case, Pair of Lead Wires, 9-Volt Battery, AC Adapter, Manual, and Pack of Four Self-Adhesive Reusable Electrodes
- HCPCS Code: E0730 for TENS Unit
- HCPCS Code: E4595 for Electrotherapy Supplies
- One-Year Warranty
- Rx Only

ORDER DETAILS:

DI8195 InTENSity™ Select Combo



TECHNICAL SPECIFICATIONS

- **Power Source:** Standard 9-volt battery or AC adapter
- **Dimensions:** 4.5" x 2.55" x 0.9"
- **Weight:** 0.28 lb
- **Channel:** Dual, isolated between channels
- **Waveforms:** monophasic square pulse, bi-phasic square pulse
- **Timer:** Adjustable from 1-60 minutes, or continuous
- **Pulse Amplitude:** Adjustable, 0-105mA and 0-70mA

STANDARD ACCESSORIES

WW3005	Lead Wires, 45"	2/pk
EP2020WC2-INTM	2" Self Adhesive Electrodes	4/pack
BA5013	9V Alkaline Battery	1 each
DI0001X	AC Adapter	1 each
CC5082	Carrying Case	1 each

