

TWINSTIM® 2ND EDITION

Portable Electrotherapy

The TwinStim® 2nd Edition is a battery-operated pulse generator that sends electrical impulses via electrodes to the body to reach the nerves and underlying muscle group. This unit is an exact combination of our InTENSity 7 and EMS 7500, which gives the user the choice of five Modes of TENS and two Modes of EMS. Features include dual channels, seven modes, timer, patient compliance meter, new stronger 100 mA output, safety amplitude cover, sturdy metal belt clip and rubber side railings for maximum handling. This device allows users a multitude of different options to address pain.

The device is provided with two controllable output channels, each independent of each other. An electrode pair can be connected to each output channel. The intensity controls are protected by a cap to avoid accidental touch. The settings are controlled by press buttons.



FEATURES:

- Offers Multitude of Pain Management Options for Nerve and Muscle Pain
- Commonly Used to Treat Neck and Back Pain
- Helps to Prevent Atrophy in Affected Muscles
- Improves Muscle Tone and Blood Circulation
- Patient Compliance Meter Can Store 60 Sets of Operation Records with a Total Recorded Time of 999 Hours
- Two Stimulation Programs: TENS, EMS
- Five Modes for TENS Therapy: Burst (B), Normal (N), Modulation (M), Strength Duration I (SD1), Strength Duration II (SD2)
- Two Modes for EMS Therapy: Synchronous, Alternate
- Intensity Level Is Controlled By the Press of a Button
- Unit Comes Complete with Device, Lead Wires, Four Self-Adhesive Reusable Electrodes, 9-Volt Battery, Hard Plastic Carrying Case and Instruction Manual
- One Year Warranty
- Rx Only

ORDER DETAILS:

DS2202 TwinStim® 2nd Edition

TECHNICAL SPECIFICATIONS

- **Channels:** Dual, isolated between channels
- **Waveform:** Asymmetrical Bi-phasic square pulse wave
- **Power Source:** One 9.0 V Alkaline battery
- **Dimensions:** 3.98" x 2.4" x 0.96"
- **Pulse Rate:** Adjustable from 2 to 150 Hz, 1 Hz/step
- **Pulse Width:** Adjustable from 50 to 300 microseconds, 10 microseconds/step
- **On Time:** Adjustable, 2-90 seconds, 1 Sec./Step
- **Off Time:** Adjustable, 0-90 seconds, 1 Sec./Step
- **Treatment Time:** Adjustable, from 1 to 60 minutes or Continuous. Adjustable in 1 minutes each step. From 1 to 15 minutes, and 5 minutes each step from 15 to 60 minutes. Treatment time countdown automatically.
- **Weight:** 150 g with battery

STANDARD ACCESSORIES

DS2202	TENS/EMS Device	1 each
WW3005	Pair of Lead Wires	1 each
EP2020WC2-INTM	2" x 2" Self Adhesive Electrodes	4/pack
TA5013-I	TENScell Alkaline Battery	1 each
	Instruction Manual	1 each
CC5082	Carrying Case	1 each
DI0009X	Wall (AC) Adapter	1 each