



# TWINSTIM® PLUS 2ND EDITION

## Portable Electrotherapy

The popular TwinStim® Plus is a four channel battery operated pulse generator that sends electrical impulses via electrodes to the body to reach the nerves and underlying muscle group. This unit is a combination simulator of TENS and EMS that can be used for pain relief and muscle stimulation. The unit comes complete with eight electrodes.

Users can now select from 4 different modes. The TwinStim® Plus gives the user the choice of five modes of TENS and three modes of EMS. The TwinStim® Plus can be powered by four AA batteries or by AC adaptor. With a stronger than ever processor, this device gives physicians and patients the capability of having a "Clinical" style device in a portable version.



### FEATURES:

- Simple to Operate, Portable Device for Pain Management
- Great Way to Prevent Muscle Atrophy
- Increase Blood Circulation
- Reduce Muscle and Nerve Pain
- Device Allows for the Use of Up to 8 Electrodes At Once
- Intensity Level Is Controlled By the Press of a Button
- Two Stimulation Programs: TENS & EMS
- Five Modes for TENS Therapy: Burst (B), Normal (N), Modulation (M), Strength Duration I (SD1), Strength Duration II (SD2)
- Three Modes for EMS Therapy: Constant, Synchronous, Alternate
- Comes Complete with Device, AC Adaptor, Four Lead Wires, Eight Self-Adhesive Reusable Electrodes, Four AA Batteries, Hard Plastic Carrying Case and Instruction Manual
- One Year Warranty
- Rx Only

### ORDER DETAILS:

**DS5002** TwinStim® Plus 2nd Edition

### TECHNICAL SPECIFICATIONS

- **Channels:** Four, isolated between channels
- **Waveform:** Asymmetrical Bi-Phasic Square Pulse
- **Power Source:** Four 1.5V AA batteries
- **Dimensions:** 5.43" x 3.07" x 1.1"
- **Pulse Rate:** Adjustable from 1 to 150 Hz, 1 Hz/step (says this in the regular technical information, but the TENS portion says 1-150Hz)
- **Pulse Width:** 50-300 microseconds, 10 microseconds/step
- **On Time:** Adjustable, 2-90 seconds, 1 Sec./step
- **Off Time:** Adjustable, 2-90 seconds, 1 Sec./step
- **Treatment Time:** Two adjustable timers, from 1-60 minutes or Continuous. Adjustable in 1 minute each step from 1 to 15 minute, and 5 minutes each step from 15 to 60 minutes. Treatment time countdown automatically
- **Weight:** 425 grams with battery

### STANDARD ACCESSORIES

|                       |                                  |        |
|-----------------------|----------------------------------|--------|
| <b>DS5002</b>         | TENS/EMS Device                  | 1 each |
| <b>WW3005</b>         | Pair of Lead Wires               | 1 each |
| <b>EP2020WC2-INTM</b> | 2" x 2" Self Adhesive Electrodes | 4/pack |
|                       | 1.5V Battery, Type AA            | 4/pack |
|                       | Instruction Manual               | 1 each |
| <b>CC5402</b>         | Carrying Case                    | 1 each |
| <b>DI0001X</b>        | Wall (AC) Adapter                | 1 each |

