

CARETEC II

Portable Electrotherapy

The Care Tec II is a two channel device that combines the pain relieving qualities of TENS with the muscle strengthening qualities of EMS therapy into one easy-to-use, versatile, and portable unit. When operating in TENS mode the unit has four modes. When operating in EMS mode the unit has three modes. An industry leader at 105 mA, this device comes standard with an AC adapter. In addition, users can easily customize placement to treat large areas (entire back) or multiple treatment sites at the same time. With the device being controlled by channel buttons, intensity levels are adjustable according to the needs of users.



FEATURES:

- Simple to Operate, Portable Device for Pain Management
- Unit Combines TENS and EMS Therapy for Optimal Treatment
- Reduce Pain and Muscle Inflammation
- Enhance the Rehabilitation of Muscles
- Increase Blood Flow and Circulation
- Two Stimulation Programs: TENS, EMS
- Four Modes for TENS Therapy: Burst, Normal, Mod 1 and Mod 2
- Three Modes for EMS Therapy: Synchronous, Asynchronous, and Delay
- Comes Complete with Carrying Case, Two Sets of Lead Wires, 9-Volt Battery, AC Adaptor, Manual and Pack of Four Self-Adhesive Reusable Electrodes
- One Year Warranty
- Rx Only

TECHNICAL SPECIFICATIONS

- **Channels:** Dual, isolated between channels
- **Waveform:** Mono-phasic square pulse wave
- **Power Source:** 9.0V DC - 1x6F22 battery or Adapter output: 9.0Vdc, 800mA
- **Dimensions:** 4.5" x 2.55" x 0.9"
- **Pulse Rate:** TENS/EMS: 1Hz-150Hz
- **Pulse Width:** 50-300 microseconds
- **On Time:** Adjustable, 1-30 seconds, 1 sec./step
- **Off Time:** Adjustable, 0-60 seconds, 1 sec./step
- **Treatment Time:** Adjustable, from 1 to 60 minutes or continuous, Adjustable in 1 minutes each step. Treatment time countdown automatically.
- **Weight:** 0.28 lb. with battery

STANDARD ACCESSORIES

DT4222	Electrical Stimulator Device	1 piece
WW3005	Pair of Lead Wires	1 each
	1.5" x 1.5" Adhesive Electrodes	4 pieces
	9V Battery	1 piece
	Instruction Manual	1 each
	Carrying Case	1 each

ORDER DETAILS:

DT4222 CareTec II