

CARETEC IV

Portable Electrotherapy

The CareTec IV Electrical Stimulator is a portable electrotherapy device featuring FOUR therapeutic modes: Transcutaneous Electrical Nerve Stimulation (TENS), Electrical Muscle Stimulation (EMS), Interferential (IF), and Russian Stim, which are used for pain relief and muscle stimulation. The stimulator sends a gentle electrical current to underlying nerves or muscles via electrodes applied to the skin. The parameters of the device are controlled by the channel buttons, with the intensity levels easily adjustable according to a user's specific needs. The portability and simple operation offers you all the options of a clinical device in a portable package.



FEATURES:

- Dual Channel Device, Allowing the Use of Four Electrodes
- Simple to Operate, Portable Device for Pain Management with Easy-to-Read Digital Display
- Strengthen Skeletal Muscles That Have Undergone Atrophy
- Reduce Pain and Relax Muscle Spasms
- Increase Local Blood Circulation
- Intensity Level Is Controlled By the Press of a Button
- Four Waveforms: TENS, EMS, IF, and Russian Stim
- Four Modes for TENS Therapy: Burst (P1), Normal (P2), Pulse Width Modulation (P3), Pulse Rate Modulation (P4)
- Three Modes for EMS Therapy: Synchronous (P1), Alternate (P2), Delay (P3)
- Three Modes for IF Therapy: P1, P2, P3, P4
- Three Modes for RUSSIAN Therapy: Constant (P1), Modulation (P2), Modulation (P3)
- Comes Complete with Carrying Case, Lead Wires, 9-Volt Battery, AC Adaptor
- Manual and Pack of Four Self-Adhesive Reusable Electrodes
- One Year Warranty
- Rx Only

TECHNICAL SPECIFICATIONS

- **Channels:** Dual, isolated between channels
- **Waveform:** Mono-phasic square pulse wave
- **Power Source:** 9.0V DC - 1x6L61 battery or Adapter output: 9.0Vdc, 800mA
- **Dimensions:** 4.5" x 2.55" x 0.9"
- **Pulse Rate:** 1-150 Hz
- **Pulse Width:** 50-300 microseconds
- **On Time:** Adjustable, 1-30 seconds, 1 sec./step
- **Off Time:** Adjustable, 1-60 seconds, 1 sec./step
- **Treatment Time:** Adjustable, from 1 to 60 minutes or continuous. Adjustable in 1 minutes each step. Treatment time countdown automatically.
- **Weight:** 0.28 lb. with Battery

STANDARD ACCESSORIES

DT4444	Electrical Stimulator Device	1 piece
WW3005	Pair of Lead Wires	1 each
	40x40mm Adhesive Electrodes	4 pieces
	9V Battery, type 6LR61	1 piece
	Instruction Manual	1 each
	Carrying Case	1 each

ORDER DETAILS:

DT4444 CareTec IV