


K3 WHEELCHAIR



K3 WHEELCHAIR

Fig.	Description	MH Part No.
K3 REPLACEMENT PARTS		
1	Probasics K3 Wheelchair Skirt Guard w/Hardware, Pair	X1WC3AS01
2	Probasics K3 Wheelchair Seat Rail Guide, Pair	X1WC3AS02
3	Probasics K3 Wheelchair Right Flip-Back Armrest Detach. Guides Front&Rear w/Hardware	X1WC3AS03
4	Probasics K3 Wheelchair Left Flip-Back Armrest Detach. Guides Front&Rear w/Hardware	X1WC3AS04
5	Probasics K3 Wheelchair 16" Seat Extender Attachment w/Seat Tube	X1WC3AS05
6	Probasics K3 Wheelchair 18" Seat Extender Attachment w Seat Tube	X1WC3AS06
7	Probasics K3 Wheelchair 20" Seat Extender Attachment w/Seat Tube	X1WC3AS07
8	ProBasics K3 Wrap Around Height Adj Arm, Pair (for SN>AS1904000547) (not pictured)	X1WC3AS08
9	ProBasics K3 Non-Wrap around Height Adj Arms, Pair	X1WCAS19
10	Probasics K3 Wheelchair Rear Wheel Axle Bolt w/Nut	X1WCAS11
11	Probasics K3 Wheelchair 16" Seat upholstery w/Hardware	X1WCAS12
12	Probasics K3 Wheelchair 18" Seat upholstery w/Hardware	X1WCAS13
13	Probasics K3 Wheelchair 20" Seat upholstery w/Hardware	X1WCAS14
14	Probasics K3 Wheelchair 16" Back upholstery w/Hardware	X1WCAS15
15	Probasics K3 Wheelchair 18" Back upholstery w/Hardware	X1WCAS16
16	Probasics K3 Wheelchair 20" Back upholstery w/Hardware	X1WCAS17
17	Probasics K3 Hand Grips, Pair	X1WCAS24
K3/K4 REPLACEMENT PARTS		
18	Probasics K3/K4 Wheelchair Right Side Wheel lock	X1WCAS03
19	Probasics K3/K4 Wheelchair Left Side Wheel lock	X1WCAS04
20	Probasics K3/K4 Wheelchair Front Caster w/Bearings	X1WCAS05
21	Probasics K3/K4 Wheelchair Front Fork Assembly	X1WCAS06
22	Probasics K3/K4 Wheelchair Stem Dust Cap - Front Fork	X1WCAS07
23	Probasics K3/K4 Wheelchair Fork Stem Bearings	X1WCAS08
24	Probasics K3/K4 Wheelchair Caster Axle Bolt w/Nut	X1WCAS09
25	Probasics K3/K4 Wheelchair Rear Wheel w/Bearings	X1WCAS10
26	Probasics K3/K4 Padded Armrest, Desk, Pair	X1WCAS20
27	Probasics K3/K4 Cross Brace, 16" - 18", Pair	X1WCAS22
28	Probasics K3/K4 Cross Brace, 20", Pair	X1WCAS23
ACCESSORY PARTS		
	Probasics K3/K4 Anti-Tippers	X1WCASANTI
	Probasics K3/K4 Elevating Leg Rest	ELR-C
	Probasics K3/K4 Swing Away Footrest	RSA-C