

PRODUCT COMPARISON

HOW TO KNOW IF THE **DELUXE** OR THE **LITE** IS THE RIGHT CHOICE?

K4 LITE



- 3-Position quick-release axles allow for tool-free removal of rear wheels for fast hemi-height adjustment
- Comes included with a seat rail and upholstery extension that easily adjusts seat depth from 16" to 18" without any additional upholstery
- Flip-back padded arm rests are height adjustable from 10" to 14"
- Adjustable height wrap-around back upholstery for maximum comfort
- Comfortable and stylish nylon upholstery
- Polyurethane tires and composite casters for superior performance
- Includes tool-free adjustable front riggings
- Adjustable anti-tippers included
- Easy push-to-lock wheel locks
- HCPCS* Code: K0195

E0990 (elevating legrests only)

E0971 (anti-tippers)

E2201 (non-standard seat frame - 20" chair only)

SPECIFICATIONS

Seat Width:	16" 18" 20"
Seat Depth:	16" – 18"
Seat-to-Floor Height ¹ :	Adult 19.5" Hemi 17.5"
Back Height:	16" – 17.5"
Back Style:	Height Adjustable
Overall Width w/Arms:	24.5" 26.5" 28.5"
Overall Depth:	32" (without riggings)
Weight Capacity:	300 lb
Product Weight ² :	32 – 36 lb
Carton Dimensions:	33" x 13" x 37"
Shipping Weight:	42 – 50 lb
Warranty:	Limited Lifetime on Frame, One Year on all Parts

vs

NEW K4 DELUXE



Scan here to
learn more about
the K4 Deluxe!

MAJOR DIFFERENCES:

- **Adjustable Angle Back** is designed to improve patient posture and positioning when using a back cushion
- **Robust lightweight aluminum frame design**
- **Convertible Armrests** are both adjustable in height and length from Desk to Full in three settings
- Stylish rear Mag Wheels have brushed aluminum handrims and quick release rear axles for fast, tool-free hemi-height adjustment and for easy breakdown and transport
- Comes included with a seat rail and upholstery extension that easily adjusts seat depth from 16" to 18" without any additional upholstery
- Adjustable height wrap-around back upholstery for maximum comfort
- The chair offers and adjustable angle caster journals that allow the rear height of the chair to lower to create "Dump" while maintaining proper front caster angle
- Adjustable anti-tippers included
- HCPCS* Code: K0004

K0195 or E0990 (elevating legrests only)

E0971 (anti-tippers)

E2201 (non-standard seat frame - 20" chair only)

SPECIFICATIONS

Seat Width:	16" 18" 20" 22"
Seat Depth:	16" – 18"
Seat-to-Floor Height ¹ :	Adult 19.5" Hemi 17.5"
Back Height:	16" – 18"
Back Style:	Height Adjustable
Overall Width w/Arms:	24.5" 26.5" 28.5"
Overall Depth:	32" (without riggings)
Weight Capacity:	16/18/20" Seat Width Chairs - 300lb 22" Seat Width Chairs - 350lb
Product Weight ² :	32 – 36 lb
Carton Dimensions:	33" x 13" x 37"
Shipping Weight:	42 – 50 lb
Warranty:	Limited Lifetime on Frame, One Year on all Wear Parts