




## Roscoe Knee Rollator

The **Roscoe Knee Rollator** offers increased comfort and ease of mobility to patients who cannot put weight on their foot or ankle, but want to remain active. With its four large wheels and padded knee platform, the Knee Rollator is perfect for indoor and outdoor use. The offset post on the aluminum knee platform allows for right-to-left adjustment in seconds - no tools required. The Knee Rollator is equipped with an easy folding mechanism that collapses using a thumb release lever and its beneath the raised knee platform for ease of storage and transportation.



Folds Flat for Easy Storage



Adjustable Hand Brake

### FEATURES

- Knee Rollator can be steered for increased maneuverability
- Ideal for individuals recovering from an ankle or foot injury
- Large 8-inch clear PU wheels with sealed bearings
- Adjustable hand brake ensures safety
- Padded knee platform with easy-adjustment mechanism
- Padded handlebar grips
- Limited Lifetime Warranty

### COLORS

Black  
Finish



White  
Finish



Burgundy  
Finish



### SPECIFICATIONS

Patient Height Range:	4'11" – 6'6"
Knee Platform Dimensions:	7" (W) x 14" (L)
Adjustable Handle Height:	31" – 41"
Knee Platform Height:	17.5" – 22"
Product Weight:	22.8 lbs.
Weight Capacity:	350 lbs.

Item #	Description	UPC	UOM
ROS-KSB	Roscoe Black Knee Rollator	092237601136	each
ROS-KSW	Roscoe White Knee Rollator	092237601471	each
ROS-KSBG	Roscoe Burgundy Knee Rollator	092237601143	each

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.