

Lumin™



## Lumin™ Multi-Purpose UV Sanitizer

FOR HOUSEHOLD ITEMS, N95 MASK DECONTAMINATION AND CPAP ACCESSORIES

The **Lumin Multi-Purpose UV Sanitizer** is the easiest and fastest way to disinfect household items, CPAP masks and accessories as well as decontaminate N95 and general use face masks without water, chemicals or harmful ozone. The Lumin's hospital-strength germicidal UV-C light source completes a disinfection cycle in just 5 minutes with a 99.9% kill rate of most bacteria and viruses, including Human Coronavirus.

UV-C light is considered to be the safest disinfection option on the market, as Ozone used in other systems has been shown by the EPA to cause irritation of the lungs and respiratory system at relatively low amounts. The Lumin Multi-Purpose UV Sanitizer is the ideal maintenance-free accessory for anyone concerned about health and hygiene.



Attractive Retail Packaging



Simple One-Button Operation



Sleek Space-Saving Profile

### FEATURES

- Originally developed out of the need for a product that would not damage the CPAP device itself during disinfection and had zero risk of exposing a respiratory patient to a lung irritant (no harmful Ozone)
- Quick 5-minute disinfection cycle time
- Hospital-grade germicidal UV light kills 99.9% of harmful bacteria, viruses, mold and fungus including Human Coronavirus
- Safe and easy to use; simple one-button operation
- Comes complete with Lumin UV Sanitizer Device, user manual and sample CPAP mask wipes
- Can also be used as a general-purpose disinfecting system for N95 and general use face masks, hearing aids, dentures, pacifiers and even toothbrushes, anything that can be placed in the drawer of the system
- Sold as each or pallet quantity of 60 units
- 24 month warranty
- Lumin is lab tested to kill over 99.9% Human Coronavirus on an N95 respirator

### PRODUCT SPECIFICATIONS

Product Dimensions:	12.25" L x 7.75" W x 8.5" H
Inner Drawer Dimensions:	9.25" L x 6.25" W x 3.75" D
Product Weight:	5.5 lbs (7 lbs in packaging)
Electrical Requirements:	110V/220V (50 - 60Hz)
Power Rating:	13 W
Operating Temperature:	20° - 100°
Operating Relative Humidity:	0 - 93%, non-condensing
Storage & Transport Temp:	0° - 120°F
MAP Price:	\$269

### PRODUCT SHIPPING SPECIFICATIONS

Shipping Dimensions:	15.25" L x 10.5" W x 10.5" H
Shipping Weight:	8 lbs

Item #	Description	UPC	UOM
LM3000	Lumin Multi-Purpose UV Sanitizer	817783020404	each