


Lumin™

Lumin™ Bullet CPAP Tubing UV Sanitizer

Introducing the **Lumin Bullet**, the newest way of sanitizing your CPAP hose, from the inside out with no chemicals! Powerful germicidal UV-C light destroys germs and pathogens in 60 seconds and stops the formation of biofilm buildup from use of heated humidification. Compatible for use on most popular CPAP tubing brands, including heated tubes, the Lumin Bullet comes complete with UV-C light, shell, USB charger and user manual.

UV-C light is considered to be one of the safest sanitization option on the market, as Ozone used in other systems has been shown by the EPA to cause irritation of the lungs and respiratory system at relatively low amounts. The Lumin Bullet is the ideal maintenance-free accessory for today's CPAP patient concerned about hygiene.



Powerful 360° UVC Light

Destroys over 99% of bacteria, mold, viruses and fungi in less than a minute. Ozone-free.



Hose Compatibility

Works with all CPAP hoses, including heated hoses. (except ResMed ClimateLine)



Retail Packaging

Attractive retail packaging includes UV-C light, Shell, USB Charger and User Manual

FEATURES

- Developed as a solution for CPAP hose sanitization and to prevent the growth of biofilm inside the wall of the hose without the risk of exposing respiratory patients to Ozone, which is known to be a lung irritant
- Quick 60 second sanitization cycle time
- Rechargeable battery lasts up to 30 uses on a single charge
- UV light kills 99.9% of harmful bacteria, viruses, mold and fungus
- Comes complete with UV-C light (Bullet), shell, USB charger and user manual
- Compatible with most popular brands of tubing in the market, including heated (not compatible with ResMed ClimateLine)
- One Year Warranty

PRODUCT SPECIFICATIONS

Product Dimensions:	8.5" L x 2.5" Diameter
Inner Diameter:	1.375" / 34.9 mm
Product Weight:	10.5 oz
Shipping Dimensions:	10.5" L x 7.5" W x 3.5" H
Shipping Weight:	2 lbs
Voltage:	DC 5V 500mA
Power Rating:	3W
Storage Temp:	0°-120°F (-18°-49°C)
Operating Temp:	20°-104°F (7°-40°C)
Operating Relative Humidity::	0-93%, non-condensing
Altitude:	Up to 6,561' (2000 M)
MAP:	\$99

Item #	Description	UPC	UOM
LM4000	Lumin™ Bullet CPAP Hose UV Sanitizer	817783020688	each