



## FEATURES

- 4 – 20 cmH<sub>2</sub>O therapy pressure
- Quiet operation with real time audible alert for tubing/mask leak
- RESlex™ exhalation relief with three comfort setting
- Remote compliance data access with optional WiFi, QR Coding, iCode™ or SD Card
- Embedded memory stores previous night's full raw data, 365 nights' user data and more than 30 years' therapy record
- SD Card data offers average statistics for specified date range with day-by-day data displaying the daily averages contained with that same date range
- Previous night's compliance-therapy snapshot to assist patient with compliance
- Maintain patient records, generate sleep reports and store patient information on iCode™ Connect, a free HIPAA-compliant patient management system
- Automatic leakage and altitude compensation
- Inspiration trigger for automatic start-up
- Two year warranty on device and humidifier

# Luna CPAP & Auto CPAP Systems

The 3B Medical Luna CPAP System with RESlex™ technology and heated humidification is a single pressure system with advanced features and is a cost effective option for treating sleep apnea. The Luna Auto CPAP System uses unique Dual-Channel Quad-Phase technology. Clinician's can set the initial and maximum allowable treatment pressures, and the Luna Auto CPAP will analyze the state of the patient's upper airway on a breath-by-breath basis and deliver pressure within the allowed range according to the degree of obstruction.

Increasing therapy comfort, RESlex™ Exhalation Relief eases the exhalation process making it easier to breathe out against the air flow. Access advanced data remotely with real-time, optional no-cost WiFi wireless data reporting, QR Coded data for smartphone upload, built-in iCode™ function or SD card, and view the data from the last session of use on the large 3.5" color LCD screen. Interactive patient coaching also available via IVR, e-mail & text. The Luna heated humidifier snaps onto the side of the device using the same power supply as the machine. Optional cellular plug-in offers the ability to capture your patient's necessary compliance data from almost anywhere, all the time.

## HCPCS CODES\*:

Device:	E0601*
Filters:	A7038*
Replacement Chamber:	A7046*
Heated Humidifier:	E0562*

## ERGONOMIC DESIGN FEATURES

- Sleek Modern Design
- Integrated heated humidification
- DC 24V powered and infrared controlled humidifier for patient safety and comfort
- Delay-off feature protects the device from humidity hazard
- Large 3.5" LCD display for easy setting adjustments
- Single dial menu access for quick setup
- Innovative dual chamber humidifier provides greater efficiency
- The most no-cost data capture options in the industry

Item #	Description	UOM
L7000	Luna CPAP System with Heated Humidification	each
L8000	Luna Auto CPAP System with Heated Humidification	each
L1020BX	Luna CPAP Air Filter	12/pk
L4510	Luna Humidifier Water Chamber	each
WL-200	Luna Cellular GPRS Module, Includes lifetime data	each
WL-100	Luna WiFi Module	each

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.