


MER-4500



REF	Item	Description
1	4502	MATTRESS
2	4502-C-88	MATTRESS COVER
3	4502V88	VENTED CELL W/T-CONNECTOR
4	PRT-4502-VCCELL	VENTED CELL W/O T-CONNECTOR
5	4502NV88	NON-VENTED CELL W/T-CONNECTOR
6	PRT-4502-NVCELL	NON-VENTED CELL W/O T-CONNECTOR
7	4581	PUMP
8	PRT-CPR-RIBBON	CPR PLUG WITH RIBBON (NOT SHOWN)
9	90165	APM MATTRESS HOSE (2 PORT)
10	PRT-CELL-CONN	MALE CONNECTOR FOR CELL (NOT SHOWN)
11	PRT-CHK-VALVE	ONE WAY CHECK VALVE (NOT SHOWN)
12	PRT-DC	CELL QUICK DISCONNECT (NOT SHOWN)
13	90121	QUICK CONNECT FITTING W/ TRANSPORT CAP (NOT SHOWN)
14	90079	7" SINGLE REPLACEMENT TUBING (NOT SHOWN)
15	PRT-4581-FILT	PUMP FILTER (NOT SHOWN)

Images are for reference purposes only and may not resemble the specific APM model.