

Ultra-Care 5800



BASIC ALTERNATING PRESSURE/LOW-AIR LOSS MATTRESS SYSTEM

The **Ultra-Care 5800** is our most economical 8" therapeutic support surface providing both alternating pressure and low air loss for the management and treatment of pressure ulcers. An ideal overlay for low risk users or the home environment, the Ultra-5800 still offers all the functionality needed from a dynamic support surface, but at a cost-effective price.



GROUP 2

- Includes mattress, cover and pump
- Pillow function maintains the air cells at the head for patient comfort
- Ventilated cell construction for low air loss
- CPR quick deflate at head of mattress
- Air cells made with nylon/PVC
- Quilted cover is water-resistant, vapor-permeable nylon; low shear/low friction
- 8 LPM pump in an easy to use design
- 10 minute cycle times
- Two modes of therapy - static and alternating
- Variable pressure setting allows for maximum comfort
- Support patients up to 350 lbs.
- HCPCS Code: E0277
- Two year warranty on mattress and pump

CELL CONSTRUCTION



- 8" alternating pressure air cells
- Air cells made with nylon/PVC
- Mattress weight: 9 lbs

4581 PUMP



- 8 liter-per-minute pump
- Front panel switch provides two modes of therapy - static or alternating
- LED backlit switches with dial on panel for simple set up and operation
- Retractable hooks for bed use

COVER



- Zipper cover is quilted vapor permeable nylon
- Water-resistant
- Easy cleaning
- Low shear/low friction

Item #	Description	Wt. Cap.	Dimensions	Qty.
MER-5800	Ultra-Care 5800 - 8" Mattress with 8 LPM Pump	350 lbs.	80" x 36" x 8"	each
5802	Ultra-Care 5800 - 8" Mattress Only	350 lbs.	80" x 36" x 8"	each

Product Specifications:

ULTRA-CARE 5800 MATTRESS

Dimensions:	80" L x 36" W x 8" H
Cells:	20 rows of cells 8" high
Inflation Time:	20 - 25 minutes with static inflate. 2 - 5 minutes with rapid inflate pump.
Material:	Nylon/PVC (do not quick-disconnect)
Weight Capacity:	350 lbs.
Mattress Weight:	9 lbs.
Cover:	Urethane coated nylon that is vapor permeable and low shear/low friction
Transport Cap:	Located at end of the connector hose which when inserted, locks and prevents air loss from the mattress if the pump loses power or is disconnected.
Additional:	CPR quick deflate at head of mattress.

ULTRA-CARE EXCEL 4581 PUMP

Input Rating:	110-120 VAC, 60Hz, 100mA, 8.5W
Power Cord:	SJT 14 Feet, 3 Prong
Air Output:	8 LPM
Indicator Lights:	Power On, Normal Pressure, Low Pressure, Mode Selection
Pressure Range:	20 - 70 mmHg \pm 10%
Cycle Time:	10 minutes
Dimensions:	11.8" L x 3.7" W x 7.5" H
Weight:	4.8 lbs.
Operating Modes:	Alternating Pressure or Static
Safety:	UL 544 & CSA

Optional Accessories:

QUICK INFLATE PUMP (ITEM #90501)



- Reduces inflation time dramatically
- 20 liter-per-minute flow pump
- Rapid deflate for quick transport
- Small and easy to use

SAFE-T-GUARD COVER (ITEM #SG8_FG)



- Perimeter cover reduces risk of falls and entrapment
- Foam borders are 4" high and sloped inward
- Zippered cover for secure fit
- Water resistant quilted top

* This is a Class II Medical Device. Federal law restricts these devices to sale by or on the order of a licensed health care professional.

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.