

Respiratory – Aerosol Therapy



SOCCER BALL NEBULIZER COMPRESSOR SYSTEM

Designed to reduce anxiety during treatments and encourage compliance, the child-friendly Roscoe Soccer Compressor Nebulizer provides effective aerosol treatments with a state-of-the-art piston compressor. Delivers 29 PSI at maximum pressure and 8 to 13 PSI operating pressure (liter flow range 5 - 8 LPM). Equipped with a convenient nebulizer parking post. Comes complete with compressor, disposable nebulizer, mask, mouthpiece, 7-foot kink-resistant tubing, five extra filters and a five year warranty. Also available with a reusable neb kit option (Model No. SOCCER-TRU & SOCCER-TRUWB-PB).



PRODUCT FEATURES:

- Child-friendly design helps to reduce anxiety during treatments
- Provides smooth, quiet treatments
- Equipped with a convenient nebulizer parking post
- Complete with compressor, nebulizer, child mask, angled mouthpiece, tubing, five extra filters and instruction manual
- TruNeb Reusable Neb Kit available for purchase with the compressor or alone
- SOCCER-TRU & SOCCER-TRUWB-PB includes Reusable Neb Kit

OPTIONAL ACCESSORIES:



Item #	Description	UOM
NEB-SOCCER	Soccer Ball Nebulizer System with Disposable Neb Kit	each
SOCCER-TRU	Soccer Ball Nebulizer System with Reusable TruNeb Kit	each
SOCCER-TRUWB-PB	Soccer Ball Nebulizer System with Disposable Neb Kit, Reusable TruNeb Kit, Carry Bag, and Piggyback Serial Labels	each
Accessories		
50016	TruNeb Reusable Neb Kit With Angle Mouthpiece And Tubing	10/cs
CIF-1128PK6	Foam Filter, 5/8 in. Diameter x 1/2 in. Thick	6/pk

PRODUCT SPECIFICATIONS:

Method of Nebulization:	Piston Compressor
Medication Capacity:	5 ml
Particle Size:	0.5 µm to 10 µm
Unit Dimensions:	5.25" L x 6" W x 7.5" H)
Product Weight:	3.3 lbs.
Maximum Pressure:	29 PSI
Operating Pressure:	8 – 13 PSI
Liter Flow:	8 LPM
Electrical Requirement:	120 VAC/60Hz
Warranty:	Limited 5 Year

HPCPS CODES:

E0570*	Nebulizer with Compressor
A7003*	Disposable Nebulizer Kits
A7005*	Reusable Nebulizer Kits
A7013*	Foam Filters
A7015*	Aerosol Mask

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

** Prescription Required: This is a Class II Medical Device. Federal law restricts these devices to sale by or on the order of a licensed health care professional.