

Respiratory - Aerosol Therapy

SIERRANEB 2 NEBULIZER SYSTEM



The Roscoe SierraNeb 2 Nebulizer System offers superior performance at an affordable price. The simplistic one-button operation and ultra-compact lightweight design makes it extremely portable and patient friendly. Weighing only 2.65 pounds, the SierraNeb 2 is easy to handle and transport, providing reliable treatment anywhere. The complete system includes compressor, nebulizer, angled mouthpiece, tubing, five extra filters, instruction manual and five year warranty. The Roscoe SierraNeb 2 Nebulizer System is also available with a reusable neb kit option.



Features:

- One-button operation
- Powerful piston pump is efficient, dependable and quiet
- Ultra compact design makes it extremely portable
- Produces a particle size of less than 5 microns
- Darker body color maintains its appearance after prolonged periods of use
- Comes complete with compressor, nebulizer, angled mouthpiece, tubing, five extra filters and instruction manual
- TruNeb Reusable Neb Kits available - purchase with the compressor or alone
- Attractive retail packaging

Specifications:

Method of Nebulization:	Piston Compressor
Medication Capacity:	6 ml
Particle Size:	Less than 0.5 μ m
Avg Nebulization Rate:	0.25 - 0.35 ml/minute
Operation Cycles:	20 minutes; off 40 minutes
Unit Dimensions:	6" (L) x 4" (W) x 3" (H)
Product Weight:	2.65 lbs.
Sound Level:	58 dBA at 1 meter
Maximum Pressure:	35 - 45 PSI
Operating Temperature:	50° - 104°
Operating Flow Rate:	12 PSI
Electrical Requirement:	120 VAC/60Hz
Power Consumption:	1.4 amps
Warranty:	Limited 5 Year

HPCPS Codes:

E0570*	Nebulizer with Compressor
A7003*	Disposable Nebulizer Kits
A7005*	Reusable Nebulizer Kits
A7013*	Foam Filters

Item #	Description	UOM
NEB-SIERRA	Roscoe SierraNeb 2 Nebulizer System with Disposable Neb Kit	each
SIERRA-TRU	Roscoe SierraNeb 2 Nebulizer System with one Disposable Neb Kit and one TruNeb Reusable Neb Kit	each
Accessories		
BAG-NEBR	Nylon Carry Bag for SierraNeb 2 Nebulizer System	each
BAG-DRW	Nylon Drawstring Carry Bag for Mini Nebulizer Compressor System, 13" x 12"	each
CIF-1128PK6	Foam Filter, 5/8 in. Diameter x 1/2 in. Thick	6/pk

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.
 ** Prescription Required: This is a Class II Medical Device. Federal law restricts these devices to sale by or on the order of a licensed health care professional.