

# truAIRE-5

by Roscoe Medical



## Stationary 5 Liter Oxygen Concentrator

Introducing the truAIRE-5 by Roscoe Medical, Compass Health Brands' first 5-liter stationary concentrator. The truAIRE-5 is a compact, low maintenance solution with easy, out-of-box set up and ultra-quiet operation. Designed to operate 24/7, the truAIRE-5 is the dependable 5-liter oxygen concentrator choice for your oxygen patients.

### PRODUCT FEATURES:

- Compact and easy to transport
- Ultra-quiet operation
- Durable, double-wobble compressor designed for 24/7 use
- Low maintenance
- Easy removal and replacement of cabinet and compressor intake filters
- 3-Year Limited Warranty

### TECHNICAL SPECIFICATIONS:

Dimensions:	13"L x 10.25"W x 21.25"H
Weight:	35.5 lbs
Liter Flow:	0.5 -5 LPM
Maximum Outlet Pressure:	5.5 psi $\pm$ 0.7 psi (35 kPa $\pm$ 5 kPa)
Oxygen Concentration:	93% ( $\pm$ 3%) at 0.5~5 LPM after warmup period of 30 minutes
Mode of Operation:	Continuous Flow
Safety System:	Low Pressure Alarm Low Oxygen Concentration Alarm High Pressure/High Temperature Shutdown Current Overload/Line Surge Shutdown
Electrical Requirements:	AC120V / 60 Hz
Rated Current Input:	3.5A
Sound Level:	$\leq$ 50.5 dBA
Warranty:	3 Year Limited Warranty

Item #	Description	UPC	Qty
O2C5L	truAIRE-5 Stationary 5 Liter Oxygen Concentrator	092237623152	each