

PBFEB3MLO Full Electric Three-Motor Bed with Standard/Low Bed Ends

NUMBER	ITEM#	DESCRIPTION
1	HBFEBEL	PB LIGHTWEIGHT FE LOW/STD HT HEAD& FOOTBOARDS W/CASTORS
	HBFEBEL	PB LIGHTWEIGHT FE LOW/STD HT HEAD& FOOTBOARDS W/CASTORS
2	X1FE3MDG06	PBFEB3M - HEAD SPRING USE WITH PBFEB3M BED
3	X1FE3MDG07	PBFEB3M - FOOT SPRING USE WITH PBFEB3M BED
4	X1HBLDG01	LOW BED HEADBOARD INSERTS (PR)
5	X1HBSMDG20	BED END PANEL INSERT - PR 90CM X 39CM, WOODGRAIN
6	X1FE3MDG01	PBFEB3M - HIGH LOW MOTOR USE WITH PBFEB3M BED
7	X1FE3MDG02	PBFEB3M - FOOT MOTOR USE WITH PBFEB3M BED
8	X1FE3MDG03	PBFEB3M - HEAD MOTOR USE WITH PBFEB3M BED
9	X1FE3MDG04	PBFEB3M - HAND PENDANT USE WITH PBFEB3M BED
10	X1FE3MDG05	PBFEB3M - JUNCTION BOX USE WITH PBFEB3M BED
11	X1HBSMDG09	HAND CRANK
12	X1HBSMDG04	DRIVE SHAFT (2 PCS)
13	X1PBHBLOCK	PBHB BED CORNER LOCKS ALL MODELS
14	X1HBSMDG03	CASTER SET (2 LOCK/2 NON-LOCK)

