




## ProBasics Satin Foam Mattress

The ProBasics Foam Mattress provides patients with the additional support and comfort they need to prevent pressure sores. Made with premium quality foam, this mattress is designed to offer both effective pressure relief and a comfortable night's sleep.

### Features

- *Premium quality foam reduces pressure on the bones to help prevent bed sores and skin ulcers*
- *Cover is waterproof and removable*
- *Non-slip underside keeps overlay firmly in place atop mattress*
- *Fire retardant and meets 16 CFR 1633 requirements*
- *Ships compressed; box is approximately 14" x 14" x 37"*
- *Weight capacity: 300 lbs/136 kg*
- *Product dimensions: 80" x 36" x 6"*
- *Limited 12-month warranty*



**Item #: PBSATINFOAM3680**

**UPC: 815067072941**

*Give patients the necessary additional comfort and support to prevent pressure sores, choose the ProBasics Foam Mattress.*