

# Precision Medical EasyPulse TOC

## Total Oxygen Concentrator

Precision Medical's Total Oxygen Concentrator is the perfect balance of mobile and stationary oxygen delivery. It delivers continuous flow and pulse flow all in the same unit, weighing less than 11.5 lbs, eliminating the need to have 2 separate oxygen concentrators. The EasyPulse TOC has an easy user interface – one knob easily sets the oxygen flow to either pulse or continuous mode. Unlike many other TOC's in the market, there are no confusing electronic controls. The convenient fold-up cart, with its adjustable telescoping handle and oversized wheels, is easily attached to the TOC. This enables the patient to take the unit anywhere and by any means. It is FAA approved for use on airlines and fits easily under an airline seat. The battery lasts over 3 hours and with optional, easy to change extra batteries, you can extend your mobility as long as you wish. All you need is in the box – included with the EasyPulse TOC is the Cart, AC Power Supply, AC Power Cord, DC Power Cord, Battery, Accessory Bag and User Manual.

### Specifications

Dimensions	7.0"L x 10.5"W x 14.63"H	
Weight	11.5 lbs (5.2 kg)	
Operating Temperature Range	TOC: 41°F - 104°F (5°C - 40°C) Battery: 32°F - 113°F (0°C - 45°C)	
Altitude	Sea Level to 10,000' (3,048 m)	
Oxygen Concentration	87% - 95% at all settings	
Sound	50.6 dBA (on Continuous Setting 2) 42.9 dBA (on Pulse Setting 1)	
Cannula	<b>Pulse Mode:</b> Maximum 7' (2.1 m) long standard adult single lumen nasal cannula <b>Continuous Mode:</b> Maximum 30' long adult standard single lumen oxygen nasal cannula	
Breath Rate Range	Up to 35 breaths per minute	
Power	<b>AC Power:</b> 100-240 VAC (1.0 A @ 120 VAC) <b>DC Power:</b> 11-16 VDC (5.0 A @ 12 VDC)	
Battery Duration (at any breath rate (approx.))	<b>Pulse Mode</b>	<b>Continuous Mode</b>
	Setting 1: 3.3 hours	1/8: 5.0 hours
	Setting 2: 3.0 hours	¼: 4.3 hours
	Setting 3: 2.8 hours	½: 3.3 hours
	Setting 4: 2.5 hours	1: 1.8 hours
	Setting 5: 2.3 hours	2: 0.8 hours
Warranty	3 Years on Device 180 Days on Battery	
HCPCS	E1390, E1392	



Item #	Description	Qty	Case QTY
PM4400	EasyPulse TOC, includes AC Power Supply, AC and DC Cords, Battery, Accessory Bag, Cart and User Manual	Each	1