



## Bariatric Shower Chair with Back

The **ProBasics Bariatric Shower Chair** is designed for individuals who require seating support when bathing or showering. The chair features a special cross brace design with a heavy duty blow-molded back and textured seat that can support patient weights up to 500 pounds. The angled legs have one-inch, easy-snap button adjustments and wide-based, non-skid feet for added sturdiness and safety. Tool-free assembly and limited lifetime warranty are standard.

### FEATURES

- Comfortable, blow-molded seat with integrated handles for additional support
- Textured contoured surface with drainage holes to reduce slipping and increase comfort
- Dual cross brace design reinforced with aircraft-style rivets for extra stability
- High-strength aluminum frame to prevent corrosion
- Legs are height adjustable in one-inch increments with heavy duty height adjustment pins to accommodate larger individuals
- Wide-based tips help increase stability
- Backrest is removable
- Quick tool-free assembly
- HCPCS Code: E0240\*

### SPECIFICATIONS

Overall Footprint:	21"W x 17.5"D
Seat Dimensions:	20"W x 12"D
Adjustable Seat Height:	16.5" – 20.5"
Weight Capacity:	500 lb
Carton Dimensions:	21" x 17" x 27"
Shipping Weight:	16.5 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	UPC	UOM
BSBCWB	ProBasics Bariatric Shower Chair with Back	815067071937	2/cs

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.