




ProBasics K2 Wheelchair

STANDARD HEMI WHEELCHAIR

The ProBasics K2 Standard Hemi Wheelchair is

designed to meet a variety of patient profiles. With a patient weight capacity of 300 pounds, a dual axle for hemi-height adjustment and multi-position caster forks, it is one of the more versatile wheelchairs on the market today. Standard features include flip back removable padded armrests, easy-to-clean vinyl upholstery, chart pocket and heel loops on the front riggings. Available with swing-away footrests or elevating legrests.







Front riggings with heel loop

wheel locks

SPECIFICATIONS

SFECIFICATIONS			
Seat Width:	16" 18" 20"		
Seat Depth:	16"		
Seat-to-Floor Height ^{1:}	Adult 19.5" Hemi 17.5"		
Back Height:	18"		
Back Style:	Fixed		
Overall Width w/Arms:	24 ½" 26 ½" 28 ½"		
Overall Depth:	31 ¼"		
Weight Capacity:	300 lb		
Product Weight ² :	41 – 43 lb		
Carton Dimensions:	33" x 13" x 37"		
Shipping Weight:	45 – 50 lb		
Warranty:	Limited lifetime on frame One year on all parts		

FEATURES

- Slide tube, powder-coated silver vein steel frame
- Dual axle provides easy transition of seat height to hemi-level
- Easy-to-clean, comfortable vinyl upholstery
- Flip-back, padded, upholstered armrests are removable
- Chart pocket on back for additional convenience
- Includes tool-free swing-away footrests with heel loops or elevating padded legrests with composite footplates
- Elevating legrests are equipped with a pneumatic piston for infinite elevating positions
- Easy push-to-lock wheel locks
- Precision sealed wheel bearings in front and rear ensure long lasting performance and reliability
- Composite handrims with 24" solid urethane tires
- 8" x 1" solid rubber front casters
- HCPCS* Code: K0002, K0001

K0195 OR E0990 (elevating legrests only)

ltem #	Seat	Arm Style	Leg Style	UPC	UOM
WC21616DS	16" x 16"	Flip Back Desk-Length	Swing-Away Footrests	815067070633	each
WC21816DS	18" x 16"	Flip Back Desk-Length	Swing-Away Footrests	815067070640	each
WC22016DS	20" x 16"	Flip Back Desk-Length	Swing-Away Footrests	815067070657	each
WC21616DE	16" x 16"	Flip Back Desk-Length	Elevating Legrests	815067070664	each
WC21816DE	18" x 16"	Flip Back Desk-Length	Elevating Legrests	815067070671	each
WC22016DE	20" x 16"	Flip Back Desk-Length	Elevating Legrests	815067070688	each

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice. NOTE: These heights can vary %-inch due to tire wear and manufacturing tolerances.

NOTE: Weights based on 18" x 16" wheelchair with desk-length arms and no front riggings. Weights will vary depending on how the wheelchair is equipped.