



Bariatric Transfer Bench

The **ProBasics Bariatric Transfer Bench** helps make transfers in and out of the bathtub safer and more comfortable. Constructed with an aluminum anodized frame, it features added cross braces for stability and a heavy duty blow-molded seat. It converts easily for left or right hand transfers. The angled legs have one-inch, easy-snap button adjustments and wide-based, non-skid feet for added sturdiness and safety. Tool-free assembly and limited lifetime warranty are standard.

FEATURES

- Aluminum frame is lightweight, durable and corrosion-proof
- Durable blow-molded bench and back have a textured surface and are contoured with drainage holes to reduce slipping
- Wide-based tips and suction cups add sturdiness and prevent slipping
- Legs have one-inch easy snap-button adjustments
- Added cross brace for stability
- Converts easily for left or right hand transfers
- Quick tool-free assembly
- HCPCS Code: E0248*

SPECIFICATIONS

Overall Footprint:	32.5"W x 21.5"D
Seat Dimensions:	25"W x 15.5"D
Adjustable Seat Height:	18" - 24.25"
Weight Capacity:	500 lb
Carton Dimensions:	33.5" x 20.5" x 17"
Shipping Weight:	28 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	UPC	UOM
BSBTB	ProBasics Bariatric Transfer Bench	815067071975	2/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.