



## Deluxe Rollator with 8-inch Wheels

The **ProBasics Deluxe Aluminum Rollator** has large 8-inch wheels for improved maneuverability and is ideal for outdoor terrain. The padded seat and removable padded backrest provide a comfortable resting area. It features easy-to-operate ergonomic hand brakes, a water-resistant nylon pouch, and a built-in fold strap for easy storage and transport. Supports patient weights up to 300 pounds.



Water-Resistant  
Nylon Pouch



Easy-to-Operate  
Ergonomic Hand Brakes





8-inch Wheels

### FEATURES

- Lightweight aluminum frame
- Two-button removable curved padded backrest
- Comfortable, flip-up padded seat
- Easy-to-operate ergonomic hand brakes
- Height-adjustable for a custom fit
- Water-resistant nylon pouch attaches under seat for added privacy and security of personal items
- Equipped with a folding mechanism for easy storage and transport
- Non-marring 8" black casters with soft grip tires, ideal for both indoor and outdoor use
- HCPCS\* Code: E0143  
E0156 (seat)

### COLORS

Burgundy Finish  Blue Finish  Black Finish 

### SPECIFICATIONS

Overall Width:	22.75"
Width b/w Hand Grips:	18.25"
Seat Dimensions:	13.75" (W) x 11.5" (D)
Seat-to-Floor Height:	21"
Weight Capacity:	300 lb
Adjustable Handle Height:	31" – 39.5"
Wheel Size:	8"
Product Weight:	15.4 lb
Carton Dimensions:	24.8" x 9.1" x 35"
Shipping Weight:	18.3 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	Color	UPC	UOM
RLA8BG	ProBasics Deluxe Rollator with 8-inch Wheels	Burgundy	815067071678	1/cs
RLA8BL	ProBasics Deluxe Rollator with 8-inch Wheels	Blue	815067071661	1/cs
RLA8BK	ProBasics Deluxe Rollator with 8-inch Wheels	Black	815067071685	1/cs

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.