

DME



## Deluxe Shower Chair with Arms

The **ProBasics Shower Chair** is designed for individuals who require seating support and stability when bathing or showering. The chair features a comfortable contoured textured seat and back that can support patient weights up to 300 pounds. Arms and back are easy to remove. The angled legs have one-inch, easy-snap button adjustments and wide-based, non-skid feet for added sturdiness and safety. Tool-free assembly and limited lifetime warranty are standard.

## **FEATURES**

- Aluminum frame is lightweight, durable and corrosion-proof
- Removable padded armrests provide comfort and added security when sitting or standing
- Blow-molded bench and back provide comfort and strength
- Quick tool-free assembly
- Wide-based crutch tips help increase stability
- Two-hole handles in seat provide security when sitting or standing
- Legs are height adjustable with one-inch snap button adjustments
- HCPCS Code: E0240\*

## **SPECIFICATIONS**

Overall Footprint:	18.5"W x 21"D	
Seat Dimensions:	16"W x 13.5"D	
Adjustable Seat Height:	17.5" – 21.5"	
Weight Capacity:	300 lb	
Carton Dimensions:	19" x 17" x 22"	
Shipping Weight:	33 lb	
Warranty:	Limited Lifetime on frame One year on all parts	

ltem #	Description	UPC	UOM
BSCWBA	ProBasics Shower Chair with Back and Arms	815067071951	4/cs

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.