



Drop-Arm Commode

The **ProBasics Drop-Arm Commode** features a durable steel frame with an easy-release push-button arm mechanism allowing for safe lateral patient transfers to and from the commode. The anti-slip, soft Hypalon® armrests provide extra comfort and support. The frame can be configured to be used as a bedside commode, toilet safety frame or raised toilet seat. The tool free removable back makes configuring the commode easy for any application.

FEATURES

- Durable steel frame provides strength and stability
- Equipped with easy-to-release arm mechanism allowing for safe lateral transfers to and from the commode
- Push-button drop-arm allows easier side transfer
- Soft, slip-resistant Hypalon® covered armrests
- Seat height readily adjusts in one-inch increments
- Bucket slides easily into place
- Removable backrest
- Comes complete with 7.5 quart commode bucket with carry handle and splash guard
- Color: white frame, gray accessories
- HCPCS Code: E0163*

SPECIFICATIONS

Width between arms:	18"
Outside width:	21"
Seat dimensions:	13.5"W x 17"D
Seat opening:	8"W x 10"D
Adjustable seat height:	17" – 22"
Weight Capacity:	300 lb
Carton Dimensions:	21.5" x 16.5" x 14"
Shipping Weight:	19.5 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	UPC	UOM
BSDAC	ProBasics Drop Arm Commode	815067071906	1/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.