



Folding Canes

ProBasics Folding Canes fold easily for storage or travel. The comfortable T-shape grip reduces hand stress and fatigue on the wrist for maximum comfort. Available in two attractive finishes with a limited lifetime warranty. Supports patient weights up to 300 pounds.



Folds compact, perfect for a purse or carry-on luggage



Contoured T-shape handle

FEATURES

- Folds in four convenient parts for easy storage
- Easy-to-use, one button height adjustment with a locking anti-rattle ring allowing for quick safe fitting
- Comfortable, contoured T-shape handle
- Easy-to-clean aluminum frame
- Adjusts in 1-inch increments
- Attractive powder-coated finish available in black or bronze
- HCPCS* Code: E0100

SPECIFICATIONS

| | |
|--------------------|---|
| Adjustment Range: | 32.5" – 36.5" |
| Cane Diameter: | 7/8" |
| Cane Tip: | 3/4" |
| Product Weight: | 13 oz |
| Weight Capacity: | 300 lb |
| Carton Dimensions: | 20.5" x 6.3" x 12.6" |
| Shipping Weight: | 16.6 lb |
| Warranty: | Limited lifetime on frame 90 days on all parts |

| Item # | Description | Color | UPC | UOM |
|--------|------------------------|--------|--------------|------|
| CNFBK | ProBasics Folding Cane | Black | 815067071241 | 6cs |
| CNFBZ | ProBasics Folding Cane | Bronze | 815067071258 | 6/cs |

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.