



Hydraulic Patient Lift

Designed and created to help make patient handling safer for caregivers by helping to prevent back injuries from manual transfers, the **ProBasics Hydraulic Patient Lift** is safe, stable and affordable. With its slim design, the ProBasics Hydraulic lift is versatile enough for use with most patient transfers. The lift features include durable steel construction with a hammertone finish, a push handle, a padded swivel bar and two locking rear casters. The ProBasics Hydraulic patient lift is light-weight, assembles and disassembles easily, and will rotate patients without side-to-side sway. Supports patient weights up to 450 pounds.



Adjustable width base with smooth low friction casters with two locking rear casters.



Rotating handle on hydraulic pump and padded bar.



Offset mast with boom style provide better lift and range.

FEATURES

- Heavy duty round tube steel construction with a chip resistant brown hammertone finish
- Six-point swivel bar attachment easily adapts to all sling styles and positions (can be used with 4 or 6 point sling)
- Low friction casters significantly improve rollability and contribute to caregiver safety. Rear two casters are locking
- 20" - 64" range allows for patient to be picked up from lying position on floor
- Pump handle can rotate from side to side for convenience of Attendant
- Offset mast and boom style provide better lift path, maximizing full travel range
- Swivel bar and pump handle are padded for comfort
- Elongated handle design serves a broad range of operator heights and offers improved handling ability
- Completely enclosed cam moves legs of adjustable base
- Base offers obstruction-free operation; ideal for cleaning
- Screw-in pin at mast base provides maximum stability and prevents removal of mast

SPECIFICATIONS

Overall Width of Base (when open):	42.25" (for maximum safety)
Overall Width of Base (when closed):	22" (to easily fit thru door openings)
Minimum Height:	22"
Maximum Height:	64"
Base Height:	6.5'
Caster Size:	5"
Weight Capacity:	450 lb
Product Weight:	68 lb
Warranty:	Limited lifetime on frame One year on all parts
HCPCS Code:	E0630

SLINGS ON NEXT PAGE ►

Item #	Description	UOM
PB0087	ProBasics Hydraulic Patient Lift	each

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



Patient Lifter Slings

MESH OR NYLON SLINGS

- Four point strap hookup
- Mesh or nylon slings are washer safe and can be used in wet or dry environments
- Commode Opening: 8.5" x 11"
- Available in a variety of sizes
- HCPCS Code: E0621



Mesh Sling



Mesh Sling with
Commode Opening



Nylon Sling

Item #	Description	Size	UOM
MESH SLINGS			
PBSL110	ProBasics Full Body Mesh Sling	Medium	each
PBSL111	ProBasics Full Body Mesh Sling	Large	each
PBSL114	ProBasics Full Body Mesh Sling with Commode Opening	Medium	each
PBSL115	ProBasics Full Body Mesh Sling with Commode Opening	Large	each
PBSL116	ProBasics Full Body Mesh Sling with Commode Opening	Extra Large	each
NYLON SLINGS			
PBSL112	ProBasics Full Body Nylon Sling	Medium	each
PBSL113	ProBasics Full Body Nylon Sling	Large	each
PBSL117	ProBasics Full Body Nylon Sling	Extra Large	each