



## Raised Toilet Seats

**ProBasics Raised Toilet Seats** are a must for anyone with bending or sitting difficulties. Available in multiple heights, with and without arms and locking options to accommodate all needs. All toilet seats are constructed of heavy duty blow-molded plastic that is durable and easy to clean.



BSRTS  
Basic Raised Seat



BSRTSL  
Locking Seat



BSRTSLA  
Locking Seat with Padded Arms

### FEATURES

- Adjustable front locking mechanism tightens easily with one hand and fits most toilets. Prevents shifting of seat when transferring (BSRTSL and BSRTSLA models)
- Elevates existing seat height by 4.5" (BSRTSL & BSRTSLA) or 5" (BSRTS)
- Extra-large hygienic hole cutout
- Tool-free assembly
- Color: white
- HCPCS Code: E0244\*

### SPECIFICATIONS

Seat Width:	17"
Seat Depth:	17"
Width Between Arms:	19.25" (BSRTSLA only)
Seat Height:	BSRTS: 5" BSRTSL: 4.5" BSRTSLA: 4.5"
Weight Capacity:	350 lb
Carton Dimensions:	BSRTS: 21" x 18" x 16" BSRTSL: 17.5" x 7.5" x 18" BSRTSLA: 21" x 8" x 18"
Shipping Weight:	BSRTS: 10 lb BSRTSL: 5.5 lb BSRTSLA: 7.3 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	UPC	UOM
BSRTS	ProBasics Basic Raised Toilet Seat	815067071968	4/cs
BSRTSL	ProBasics Raised Toilet Seat with Lock	815067071982	1/cs
BSRTSLA	ProBasics Raised Toilet Seat with Lock and Padded Arms	815067071975	1/cs

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.